


Using the Rolyan® Heel Boot

(NOTE: Heel Boot is packaged the way it's applied, not inside-out)

- 1** Place heel in hole (NOTE: it may be helpful to bend the patient's knee slightly).



- 2** If also using a compression device, route tubing through port hole, ensuring that tubing neither kinks nor is in contact with patient's skin.



- 3** Tighten black straps for a snug fit, but do not over-tighten such that blood flow is restricted.



- 4** Place foot in neutral or slight dorsiflexion and secure blue foot positioning straps.



- 5** Feel underneath to make sure the patient's heel is floated. If not, undo straps and re-apply.



- Use only on non-ambulating patients. Do not allow patient to walk or stand while wearing Rolyan® Heel Boot
- These splints are to be fitted initially by a healthcare professional who is familiar with the purpose for which they are prescribed.
- The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.
- If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.
- Be careful that the straps are not so tight that they interfere with circulation.
- Recommended for use by only one patient.