## Rolyan® Preformed Humerus Brace:


A511-LS through LL and A511-RS through RL

## Rolyan® Preformed Ulna Brace:

A623-1 through 3

A

B

Be sure to note which brace you are using and follow the instructions that apply to that brace. Unless otherwise stated, these instructions apply to both of the braces.

### INDICATIONS

- To stabilize fractures

### CONTRAINDICATIONS

- Not for use over open wounds and infected areas

### PRECAUTIONS

- **These braces are to be fitted initially by a healthcare professional who is familiar with the purpose for which they are prescribed. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.**
- **If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**
- **Brace should fit snug enough to prevent migration, but not so tight that circulation is impaired.**
- **After heating splinting material, always check its temperature before applying it to the patient.**

### CARE AND CLEANING

- The brace will lose its shape in temperatures over 135°F (57°C). It should be kept away from sources of heat such as ovens, hot water, open flames, and sunny car windows.
- The brace can be cleaned with soap and lukewarm water. The straps can be washed with soap and water, but the adhesive portion should not be moistened. The brace and straps should be allowed to dry thoroughly before reapplication of the brace.

### PROCEDURE

1. If desired, cover the affected area with stockinette.

#### Preformed Humerus Brace

The proximal end has a continuous curve; the distal end is shaped irregularly.

2. Open the straps. Position the brace onto the arm so the opening and the polyethylene tongue that covers it are anterior and the brace is high into the axilla.
3. Slide the tongue under the thermoplastic shell and fasten D-ring straps. A properly applied brace is shown in **A**.
4. If modifications in the fit of the brace are necessary, see instructions 10 and 11.

#### Preformed Ulna Brace

The narrower end should be placed distally.

5. Apply the self-adhesive straps as appropriate.
6. Position the brace with the opening on the volar surface of the forearm.

7. Slide the padded polyethylene tongue under the two edges of the brace.

8. Secure straps. A properly applied brace is shown in **B**.

9. If modifications in the fit of the brace are necessary, see instructions 10 and 11. The brace can be modified to prevent pronation and supination as follows:

—Remove straps

—Heat dorsal portion of brace (see instruction 10).

—Put brace onto forearm and press a groove about ½" (1.3cm) wide over the full length of the brace to separate the interosseous membrane.

#### Modification in the fit of both braces

10. These braces are made from ¼" (3.2mm) **Polyflex II\*** splinting material with a **Kushionflex\*** lining. To trim or modify the brace, spot heat it with a heat gun or in a heat pan of water at a temperature of 150° to 160°F (65°–70°C) for one minute or until pliable. **Heating the entire brace may cause it to lose its preformed shape. Do not heat where straps are attached.** For better results when modifying the material, shape it with stroking movements.
11. The brace can be trimmed with scissors. The edges should be flared and smoothed.

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**Patterson Medical Ltd.**