


Impulse[®] 3000T TENS



BioMEDICAL
LIFE • SYSTEMS

Impulse[®] 3000T



SPECIFICATIONS

Size	3.8" x 2.60" x 1" (97 mm x 62 mm x 25 mm)
Weight	3.1 oz (87 g)
Power Source	9V Battery, E-block, type 6F22
Channels	Dual
Waveform	Asymmetrical, biphasic square-wave
Amplitude	0 – 98 (mA)
Frequency	2 - 150 (Hz or pps)
Pulse Width	50 - 250 microseconds (µs)
Output	Constant current
Continuous	Continuous stimulation
Burst	8 pulses per burst/2 bursts per second
Modulation	Decrease/increase of set width of 50% over a five-second period
Timer	15, 30, 60 minutes
Number of Electrodes/ Leadwires	4
Tolerances	+/- 10%

Data recorded across a 500 OHM load resistance

SAFETY

Caution Federal law (USA) restricts this device to sale by or on the order of a physician so licensed by the State.

Indications Transcutaneous Electrical Nerve Stimulation (TENS) devices are used for the symptomatic relief and management of chronic (long-term) intractable pain and as an adjunctive treatment in the management of post-surgical and post-traumatic acute pain problems.

Contraindications TENS devices can adversely affect the operation of demand-type cardiac pacemakers. TENS is not recommended for patients with known heart disease without a physician's evaluation of risk. Do not stimulate over the eyes or carotid sinus nerves. Do not apply TENS for undiagnosed pain syndromes until etiology is established. Do not place electrodes in a manner that causes current to flow transcranially (through the head).

Warnings This device should be used only under the continued supervision of a physician, or outside the USA, by a qualified pain management specialist. TENS is ineffective for pain of central origin, (i.e. appendicitis, hepatitis). TENS is of no curative value; it is a symptomatic treatment which suppresses pain sensation which would otherwise serve as a protective mechanism on the outcome of the clinical process. Safety of TENS devices for use during pregnancy or delivery has not been established.

Electronic equipment such as ECG monitors and ECG alarms may not operate properly when TENS is in use. The user must keep the device out of the reach of children. TENS is for external use only.

Precautions Avoid adjusting controls while operating machinery or vehicles. Turn the stimulator off before applying or removing electrodes. Isolated cases of skin irritation may occur at the site of electrode placement following long-term application. Use only for the specific pain problem as prescribed by the physician, or outside the USA, by a qualified pain management specialist. Effectiveness is dependent upon patient selection by a qualified pain specialist.

Adverse Reactions Possible allergic reaction to tape or gel. Possible skin irritation or electrode burn.

EQUIPMENT not suitable for use in the presence of a FLAMMABLE ANESTHETIC MIXTURE WITH AIR or WITH OXYGEN OR NITROUS OXIDE.

BioMEDICAL[®]
LIFE • SYSTEMS

P.O. Box 1360
Vista, California 92085-1360
(1) 800-726-8367
FAX: (1) 760-727-4220
www.bmls.com
information@bmls.com

DISTRIBUTED BY: