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Deep-Heat™ Therapy:

Relieve serious joint and muscle pain for hours with professional-quality intense moist heat!

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Home Heat Therapy:

Effective moist heat therapy with a value price-point!

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- ☒ Large 13" x 13"
- ☒ Joint Wrap 5" x 16"

Gentle Overnight Heat Therapy:

Minimize the discomfort of cold, painful feet and legs!

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- ☒ BedWarmer™
- ☒ Deluxe BedWarmer™

Light Exercisers:

Convenient, portable exercise for virtually anyone!

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2878 — 4/2011
Ice It!® HTUB - FR

Enhancing Comfort, Enhancing Treatment, Enhancing Lives!™

Instructions

Available in versatile Multi-Use styles, custom-shaped Joint-Specific styles, as well as Refill Kits and individual cold packs.

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ICE IT!®
ColdCOMFORT™ System

Targeted, convenient cold therapy to get you back in the game of life!

307 W. Jackson St. / Battle Creek, MI 49037-2369
800-253-0854 www.BattleCreekEquipment.com

Ice It![®] ColdCOMFORT™ products help you be prepared for life's little emergencies at a moment's notice. Just store the flexible cold pack in your freezer so it's always ready to use, since immediate treatment with cold therapy can greatly reduce pain and trauma after injury.

Many people overlook the effectiveness of safe, simple cold therapy to reduce tissue damage, swelling and pain caused by sprains, bumps, bites or bruises. Cold therapy can also alleviate the pain of arthritis, relax muscle spasm and help fight inflammation by increasing circulation, and is often prescribed after medical or dental procedures. Used alone, or to enhance the benefits of medication, exercise or massage, cold therapy is a powerful pain management tool.

Application: Simply remove as many frozen Ice It![®] packs from the freezer as needed and insert directly into the ColdCOMFORT™ holster*. Position over affected area and fasten in place by adjusting the Velcro closures of our Comfort-Flex™ design for gentle compression. The frozen Ice It![®] pack stays soft and pliable, to mold around painful areas with no need for you to hold anything in place.

*If you do not have our soft, insulated ColdCOMFORT™ system, do NOT place the Ice It![®] pack directly against your skin.

Application of cold packs should not exceed 20 minutes at a time. Allow 20 minutes between applications. Professionals also recommend elevating swollen areas (ankle, knee, elbow, etc) if possible, to enhance the benefits of cold therapy. Use cold packs during the first 48 hours after injury to reduce pain and swelling, or as recommended by your physician.

Caution: Persons with circulatory problems should consult their physician before using any cold pack. Limit applications to 20-minute intervals to avoid skin damage. If pain and swelling persist or increase, please consult your physician.

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