HB-RGB

INSTRUCTIONS

HAND BASED - RADIAL GUTTER BRACE

antimicrobial treated

HAND BASED RADIAL GUTTER BRACE

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



no. 10097 Rev D
USA and International Patents Pending

Medical Device & QA Services Ltd
76, Stockport Road, Timperley, Cheshire, WA15 7SN,
United Kingdom
Tel +44 161 870 6751 Fax +44 161 903 9787
Email: info@mdqservices.com
www.mdqservices.com

ENGLISH

Before using the device, please read the entire instruction manual. Correct application and care is vital to the proper functioning of the device.

INTENDED USE / INDICATIONS

The Hand Based Radial Gutter Brace (HB-RGB) provides stabilization to the second and third metacarpals and metacarpophalangeal (MCP) joints, and the second and third digits.

CONTRAINDICATIONS

Patients with known allergies to foams, adhesives, or plastic products.

⚠️ WARNINGS and PRECAUTIONS

- ▶ Brace should be fitted by a qualified healthcare professional.
- ▶ Do not heat brace on the body.
- ▶ Check to make sure brace is not too hot prior to application.
- ▶ Brace should be appropriately sized to the patient to ensure proper support and function.
- ▶ Care should be taken with patients with delicate or compromised skin, especially the elderly, not to apply the brace too tightly, which MAY cause pressure and shear to the fragile tissue. PATIENTS SHOULD BE INSTRUCTED ON APPROPRIATE TENSIONING DURING ONGOING USE. An under sleeve or STOCKINETTE may be used for added protection.
- ▶ Do not heat Exos™ products to greater than 212°F (100°C).
- ▶ FOR SINGLE PATIENT USE ONLY.

HEATING AND APPLICATION

- ▶ Proper heating temperature is 200°F (93°C).
- ▶ Heat brace in Exos oven approximately 5 minutes or until brace is pliable. Brace should not be too hot to handle comfortably.
- ▶ Ensure that skin is completely dry before application.
- ▶ Align brace on palm so curve follows the Thenar crease. (Figure 1)
- ▶ Wrap brace around hand. (Figure 2)
- ▶ Secure hook-and-loop fastener straps as shown. Two straps may be used to secure the fingers, and two straps secure the wrist. (Figure 3)

Note: Larger braces are packaged with slightly longer straps intended to be used to secure the wrist.

▶ Mold brace into desired position. (Figure 4)

Warning: The brace should not be worn too tightly. This may produce excessive mechanical force or shear forces to the skin.

CLEANING INSTRUCTIONS


1. Loosen or remove the brace.
2. Wash inside and outside of brace and skin with antibacterial soap and water (86°F / 30°C).
3. Rinse thoroughly with water.
4. Use a hair dryer set on high volume and cool setting to thoroughly dry the interior of the brace, and skin.
5. Retighten the brace.

For fitting tips on each Exos™ device: Visit the Exos™ Provider Center at: <http://exosmedical.com/provider/HB-RGB> or Scan the QR image with your smartphone

MATERIAL USED
Shell: EVA foam, spandex, nylon, terpolymer, polyester
Antimicrobial treatment

LIMITED PRODUCT WARRANTY
Exos™ covers defects in workmanship and materials up to 8 weeks from date this product is dispensed. Damage caused by normal wear and tear, abuse, or alterations is not covered.

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