

Snug Hug



Unfold blanket and set out in seating area. Separate the straps and lay them open.



1. Sit the user in the seat. Wrap the Leg Hug around the thighs.



2. Wrap the Hip Hug around the hips.



3. Wrap the Body Hug on the torso.



4. Wrap the Shoulder Hug straps around each shoulder.



5. Re-tighten all straps as needed. Wrap the blanket around the user for discretionary comfort.

Care Instructions

Snug Hug: Hand wash in cold water with mild detergent. Lay flat to dry.

Fleece Blanket: Cold water wash, delicate cycle. Air dry.

More Information

Adjusting to the user, Folding for on the go, taking apart for Washing

<http://tinyurl.com/UsingSnugHug>