

Days™ Lightweight Rollator Walker

0815-61-661

days™
a Patterson Brand



The Days Lightweight Rollator Walker is a multipurpose walker manufactured from painted steel tubing with four solid tyre wheels. It incorporates a rest seat, basket and locking brakes. This walker also has adjustable handles and can be folded flat for storage or transportation.

Maximum user weight: 364 lbs.
Handle height adjustment: 34.25"–39"

Assembly Instructions

- 1) Remove all packaging.
- 2) Fit cross bar into front sockets on the frame and push fully down.

To secure the cross bar make sure the retaining spring buttons have gone all the way through the hole in both sockets.

Handle Adjustment

Remove the hand knob and bolt. Pull the handle to new setting (there are five fixed settings at 1" intervals), re-fit and tighten the knob and bolt ensuring that the head of the bolt is inside the hexagonal hole.

Operation

The Walker is fitted with locking brakes. To apply the brakes, pull up on the brake handles. To lock the brakes, push down on both brake handles until they "click" into the locked position. To release the brakes, simply pull up on the brake handles until they "click."

Before using the walker as a rest seat, both brakes must be applied and the frame pushed forward to check that it is not moving. The brakes only lock the back wheels so the user must not expect the frame to be completely immobile.

If the user finds the brakes hard to apply, the pressure can be adjusted by reducing the tension in the spring above the brake shoe. Hold the hexagonal nut and turn the knurled nut to vary the brake pressure i.e. the distance between the shoe and the tyre, when the brakes are off. It is a case of trial and error to suit individual users' grip strength.

WARNING! Brakes that are adjusted to give light pressure will be less effective at preventing rear wheel movement under load.

When using the walker, be aware that stability depends on all four wheels being in contact with the flat surface. Care should be taken when negotiating uneven ground, inclines, curbs, etc. as stability will inevitably be reduced.

Folding for Storage or Transportation

To fold the walker for storage or transportation pull on the strap, thus pulling the locking mechanism up and allowing the wheels to come together. Pushing down on the handles unfolds the walker.

Precaution

Daily, before use, check the effective operation of the brakes. Check that the wheels turn and swivel freely and there is no debris building around the axle. Periodically check that all screws, nuts and bolts are securely fastened. Check moving parts for wear.

Cleaning and Maintenance

Minimum maintenance: occasionally wipe over with a soft cloth. Do not use harsh abrasive materials.

Limited Product Warranty

Patterson Medical warrants that all products shall be free from manufacturing defects for a period of one year from the date of shipment. OTHER THAN THE FOREGOING WARRANTY, PATTERSON MEDICAL MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, AS TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER MATTER RELATING TO THE PRODUCT. Customer must notify Patterson Medical of any manufacturing defect within five days of discovery of the defect. Upon a determination by Patterson Medical that the product is defective, Patterson Medical sole and exclusive obligation shall be to provide customer with a replacement product or a full refund of the purchase price, at the sole discretion of Patterson Medical, and to reimburse customer for any return freight costs. In no event shall Patterson Medical be liable for any other costs or expenses including lost profits or direct, indirect, consequential, incidental or punitive damages. The warranty set forth herein shall govern all purchases, notwithstanding any inconsistent terms and conditions of customer's purchase order.


PATTERSON
MEDICAL
www.pattersonmedical.com

© May 2011 Patterson Medical Products, Inc.
081561661LIT

Patterson Medical
1000 Remington Blvd., Suite 210
Bolingbrook, IL 60440
Phone: 800-323-5547
Direct Phone: 630-378-6000
Fax: 800-547-4333
Direct Fax: 630-378-6310