

## HOME REHABILITATION PROGRAM

### Neuromuscular Electrical Stimulation (NMES) Chattanooga® Continuum™

1. Clean skin thoroughly with soap and water.
2. Connect lead wires.
3. Place electrodes on skin as instructed by your health care professional.
4. Turn machine on by pressing the grey Power button in the bottom middle of the device.
5. Set device to:



Program \_\_\_\_\_

(NMES - Muscle reeducation/strengthening) \_\_\_\_\_ times per day \_\_\_\_\_ min per session

Program \_\_\_\_\_

(TENS - Pain) \_\_\_\_\_ times per day \_\_\_\_\_ min per session

6. Increase intensity by pressing upper buttons on the left and right side of the screen to reach desired level.
7. Turn off by pressing the grey power button on the bottom middle of the device.



"Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Information related to various health, medical, and fitness conditions and their treatment is not meant to be a substitute for the advice provided by a physician or other medical professional."  
 Copyright © [2019] by DJO, LLC All rights reserved. No part of this publication may be reproduced, distributed, or transmitted in any form or by any means.



**DJO, LLC** | 1430 Decision Street | Vista  
 CA 92081-8553 | U.S.A.  
[www.DJGlobal.com](http://www.DJGlobal.com)