Symbols glossary located at:

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www.performancehealth.com/medical-device-glossary

Distributed by:

Performance Health Supply Inc.

28100 Torch Parkway, Ste. 700
Warrenville, IL 60555, USA
Phone: 800-323-5547
Fax: 800-547-4333

Distributed by:

Performance Health Canada, Inc.

6675 Millcreek Drive, Unit 3
Mississauga, ON L5N 5M4
Phone: 800-665-9200
Fax: 877-566-8611

Manufactured for:

Performance Health Supply Inc.

W68 N158 Evergreen Blvd.
Cedarburg, WI 53012, USA
Phone: 800-228-3693

CEpartner4U
Esdoornlaan 13
3951 DB Maarn
The Netherlands

Country of Origin: Taiwan
Rev. 3, Effective: April 2019
©Performance Health
7101865

Reach Massager

7101865

The Rolyan Reach Massager is designed to allow users to treat muscle dysfunction and achieve effective pain management by applying direct pressure without assistance.

- Helps users independently reduce painful knots in muscles
- Makes it easy to apply deep pressure massage without assistance, even to hard-to-reach muscle groups
- Facilitates stress and tension relief without drugs via precision control
- Provides soothing, penetrating massage for faster muscle recovery
- Ergonomic design features 8 strategically-positioned treatment balls and facilitates optimum leverage
- Compact, durable & lightweight for easy transport

INITIAL USE

When first using the Rolyan Reach Massager, muscles will require a “breaking in” period, so begin with light pressure just a minute or two a few times per day.

It is important to apply adequate pressure to keep balls in contact with the muscle being worked on. Applying adequate pressure also prevents massager from sliding.

Below are different methods of applying pressure to help you determine which is most comfortable and effective for the specific muscles you are massaging:

- For extremely sensitive tender points, massage the surrounding area first.
- Thicker, deeper muscles (such as back muscles) require more pressure.
- Apply steady pressure for 5 or 10 seconds, release, relax and repeat.
- Maintain direct pressure while moving the ball slowly over a small area of the muscle, gradually going in deeper.
- Gently wiggle the Rolyan Reach Massager to burrow into a sore area.

Should you apply too much, or too prolonged pressure, you may experience discomfort. As you become more familiar with, and accustomed to working with the Rolyan Reach Massager, you will be able to massage with as much pressure as you desire and for extended periods without experiencing discomfort.

CAUTION

Do not attempt to use the The Rolyan Reach Massager for treatment of serious medical ailments, or injuries. Rolyan and its affiliates will not be held liable for misuse of the Rolyan Reach Massager in attempted treatment of unknown origins. Please consult with a qualified medical practitioner for treatment of any serious medical condition.

- When first using the Rolyan Reach Massager, limit sessions to short periods. Overuse will likely result in soreness for a day or two.
- Do not press directly on or against the spinal column or sciatic nerve.
- Avoid massaging areas of inflammation.
- If bruising occurs, discontinue use.
- Keep out of reach of children. Do not allow them to attempt to use or play with the Rolyan Reach Massager.

CARE Wash with mild soap and water.

LIMITED PRODUCT WARRANTY

Performance Health warrants that all products shall be free from manufacturing defects for a period of one year from the date of shipment. OTHER THAN THE FOREGOING WARRANTY, PERFORMANCE HEALTH MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, AS TO MERCHANTABILITY, NONINFRINGEMENT, FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER MATTER RELATING TO THE PRODUCT. Customer must notify Performance Health of any manufacturing defect within five (5) days of discovery of the defect. Upon a determination by Performance Health that the product is defective, Performance Health's sole and exclusive obligation shall be to provide customer with a replacement product or a full refund of the purchase price, at the sole discretion of Performance Health, and to reimburse customer for any return freight costs. In no event shall Performance Health be liable for any other costs or expenses including lost profits or direct, indirect, consequential, incidental or punitive damages. In the event the product includes a warranty card or extended warranty from the manufacturer, customer's sole recourse under such warranty shall be against the manufacturer. The warranty set forth herein shall govern all purchases, notwithstanding any inconsistent terms and conditions of customer's purchase order.

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