



Marena Reusable Antibacterial Face Mask*

Style: FMRAB

Washable, reusable fabric face mask with Silvadur™ antibacterial protection to guard face, nose and mouth against respiratory droplets and splashes.

This mask provides a physical barrier protection against larger respiratory droplets, germ particles, and splashes for the user and the environment. While this mask will provide some protection against, it will not act as a barrier against smaller respiratory droplets or germ particles.

HOW TO PUT ON:

When putting mask on, place shiny side to face, with pleats facing downward. Put loops over ears, and adjust to tightly cover the nose and mouth (or respirator if combined). If a tighter fit is required, the loops can be twisted before placing over ears.

HOW TO SAFELY REMOVE:

To remove contaminated mask, avoid touching outer layer of mask or your face. Remove from ears, only touching loops, and immediately put into washer, sterilizer or safe storage.

HOW TO CLEAN:

Clean mask by washing or sterilizing before re-use every time. Sanitize hands immediately after handling contaminated masks each time. Clean after multiple use or if soiled. Discard if damaged or torn.



Machine wash; tumble dry

Do not use fabric softener; Do not use bleach

This face mask is washable multiple times without losing effectiveness, even at high temperatures – $\geq 160^{\circ}\text{F}$ ($\geq 71^{\circ}\text{C}$).

May be steam cleaned at temperatures between 212°F (100°C) and 270°F (132°C).

Ironing is permitted at temperatures up to 365°F (185°C).

*This product is non-sterile and has not been tested to N95 standard. It is not a medical device. It is not intended to be used as a respirator. Each mask is intended for use by one person only and must not be shared. This mask has not been tested for effectiveness against the COVID-19 virus.

For more information, contact us at info@marenagroup.com or call 1-888-462-7362.