


Rolyan* Palm Protector with Finger Separators

A812-01 and A812-02




Strap fastens to hook



Strap sewn to splint



B



C



D

INDICATIONS

- May be used for moderate to severe finger spasticity or contractures.

CONTRAINDICATIONS

- Not for use if uncontrolled edema or open sores are present.

PRECAUTIONS

- The splint should fit snugly and comfortably without interfering with circulation.
- If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued.

CARE & CLEANING

- Fasten straps. Hand wash fabric in warm water and mild soap. Do not wring. Dry at room temperature. Be sure splint is completely dry before reapplying it.
- Unthreading the finger separating strap from the splint before washing will reduce drying time. To rethread the strap, see A.

INSTRUCTIONS FOR USE

1. With the straps open, place the fingers into the finger separators (B).
2. Place the thumb through the thumb hole (C).
3. Attach the finger strap around the palmar surface of the fifth digit (C).
4. Attach strap on dorsum of wrist.
5. Adjust tension of straps for a proper fit.
6. A properly applied splint is shown in D.
(Design copyright Smith & Nephew.)