






-  **Frequency** – Two or more times monthly. Strava recommends each facility determine its cleaning frequency based on usage, water grade, visual inspection, and sanitation protocols that apply to your facility and/or JCAHO.


 **1 TURN UNIT OFF FROM THE SWITCH ON THE BACK OF THE UNIT AND EVACUATE WATER**  
See [User Manual](#) procedure.


 **2 ADD WATER & CLEANING SOLUTION**  
Fill unit just above the water sensor with a 2:1 ratio of tap water to cleaning solution.


 **Cleaning Solution = Vinegar -OR- Coffee Machine Commercial Descaler**

 **3 TURN THE UNIT BACK ON** and leave it running for 20-30 minutes to allow the water/cleaning solution mixture to heat.


 **4 EVACUATE WATER**  
Wait 10 minutes then evacuate water. See [User Manual](#) procedure.

 **5 TURN UNIT OFF**  
Unplug and allow unit cool down before proceeding.

 **6 CLEAN BI-FOLD LID & INSIDE WATER BASIN**  
Using a Scotchbrite® Scour Pad or cloth, wipe the inside tub with warm water. To remove scum, waterlines, mineral deposits, etc., add a mild detergent. Be sure to rinse thoroughly and wipe dry. **Do NOT use bleach, ammonia-based cleaning products, or extra abrasive products such as steel wool.**

 **7 REFILL WITH WATER**  
Fill unit to operating level with tap water. See [flip side](#) or [User Manual](#) if you use distilled water.

 **8 TURN UNIT ON**

 Repeat procedures if necessary

**See Flip Side For Annual  
Maintenance Recommendations**



## SPLINT BATH CARE CARD

# QUICK TIPS

### DISTILLED WATER

Tap water is required for the low water sensor to work properly.

Follow these proportions if you are in a hard-water environment and require distilled water.

1 Cup Tap Water



- PER -



1 Gal.  
Distilled

### WATER EVACUATION

- Turn unit off using the power switch at the back.
- Allow the water to cool before water evacuation.
- Be sure that the supplied drain hose is properly inserted at the back the unit.
- Turn the quarter turn valve counter-clockwise to release the water.



### ABOUT US

Our team at Rolyan consists of seasoned executives with extensive backgrounds in the fields of Rehab, Physiotherapy, and Occupational Therapy. With years of experience under our belts, we are dedicated to introducing groundbreaking products to these markets, aimed at improving the experiences of both therapists and patients. Our primary goal is to continually innovate and refine our product lineup to meet the ever-changing demands of healthcare professionals. We take great pride in the outstanding quality of our offerings, all of which are crafted with care right here in the USA.



### FOR TECHNICAL SERVICE

Visit [www.stravasolutions.com](http://www.stravasolutions.com) and submit a Service Inquiry



Refer to the user manual at: [www.stravasolutions.com](http://www.stravasolutions.com)

PART NUMBER 400.613 ROLYAN SPLINT BATH CARE CARD

©Copyright 2024 Pivotal Health Solutions All rights reserved