(EN) **GLUTE BAND LOOPS**

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**Enable high impact glute
activation and lower body fitness**

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Four increasing resistance levels

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Thank you for purchasing
THERABAND® Glute Band Loops,
the versatile, easy-to-use exercise loops that are
ideal for use at home or while traveling.

Please visit **THERABAND.COM** for the most up
to date instructions and information.

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- Before beginning any type of exercise program, read all exercise descriptions included with the product.
- Consult your healthcare provider before starting any exercise program. If you have a medical condition, seek their approval before using.
- For beginners, perform exercises without Band Loops until you are comfortable, and then begin exercises with the Band Loops.
- Begin with one set of 10-15 repetitions of each exercise. Rest between sets. Gradually increase your sets and repetitions.
- When you are comfortable performing three sets of 10-15 repetitions, gradually progress to the next color of THERABAND® Glute Band Loop.

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EN

MONSTER WALK + SIDE STEPS

1. Loop band just above your knees.
2. Take an athletic posture, knees slightly bent, tight abs, with head and feet facing forward.
3. Point feet forward, slightly wider than hip distance apart.
4. Walk forward, taking short steps, keeping your knees apart for 10 steps.
5. Facing the same direction, move back 10 steps to the starting position.
 - Also perform **Side Steps**: move one foot directly to the side, leading with your knee, not your ankle.
 - Follow with the other foot and repeat after 10 steps.

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EN

CLAMSHELL/ SIDE-LYING

1. Loop band just above your knees.
2. Lay on your side with knees bent at 90 degrees and hips bent to 45 degrees.
3. Lift your top knee while pushing your bottom knee down to floor.
4. Keep ankles together. Squeeze your glutes. Hold and slowly return.

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(EN)

STANDING KICK BACK + SIDE KICK

1. Stand with feet hip-width apart. Place band just above knees and use wall to maintain balance.
2. Controlling movement from the glutes, kick one leg back, keeping it straight at all times.
3. Squeeze your glutes at the top of movement and slowly return to starting position.
 - Also perform **Side Kicks**: Abduct one leg to the side keeping it straight at all times.
 - Squeeze your glutes at the top of the movement and slowly return to starting position.

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(EN)

DONKEY KICK

1. Loop band just above your knees.
2. Go to your hands and knees with your back neutral and abs braced.
3. Keep your face down and neck in line with your spine, and do not arch your back.
4. Drive one leg up while maintaining 90-degree angle. The sole of your foot should face the ceiling.
5. Lower slowly to starting position.

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(EN)
FIRE HYDRANT

1. Loop band just above your knees.
2. Go to your hands and knees with your back neutral and abs braced.
3. Keep your face down and neck in line with your spine, and do not arch your back.
4. Drive your knee out to the side maintaining 90-degree angle.
5. Lower slowly to starting position.

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(EN)
GLUTE BRIDGE

1. Place band just above your knees.
2. Lie on back with knees bent, feet hip-width apart.
3. Brace abs.
4. Thrust hips up while driving knees out. Squeeze glutes and go as high as possible. Hold for four seconds.
5. Slowly lower hips to ground.

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(EN)
**SIDE-LYING
ABDUCTION**

1. Place band just above your knees.
2. Lie on your side with lower knee bent. Keep your spine straight and abs braced.
3. Straighten upper leg and turn foot inward.
4. Raise upper leg as high as possible and slowly return to starting position.

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(EN)
**90-DEGREE
CLAMSHELL**

1. Loop band just above your knees.
2. Lie on your side, bending your hips to 45 degrees and knees to 90 degrees.
3. Ensure hips are straight, knees are stacked, and feet are aligned with your back.
4. Lift upper leg as high as possible. Squeeze your glutes.
5. Return to bottom position. Your feet and knees should not touch between reps.

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
The Hygenic Company LLC

1245 Home Avenue
Akron, OH 44310 USA
(330) 633-8460
www.THERABAND.com

Distributed in the United Kingdom By: Performance Health International LTD

Nunn Brook Road
Huthwaite, Sutton-in-Ashfield
Nottinghamshire, NG17 2HU, UK
+44 (0) 1623 448 706
www.performancehealth.co.uk

Distributed in Canada By: Performance Health Canada, Inc.

2810 Argenta Road, Unit 2
Mississauga, ON L5N 8L2
Canada
+1800-665-9200
www.performancehealth.ca

**Please visit THERABAND.COM for the most up to
date instructions and information.**

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Always read the label. Follow the instructions for use. THERABAND® the Color Pyramid and Associated Colors are all trademarks of Performance Health and/or its subsidiaries and may be registered in the United States and other countries. © 2024 All rights reserved. Unauthorized use is strictly prohibited. P11978-R00