

Stall Bars

Directions

Begin by placing the rungs into the oval areas of the uprights. Use upright 'A' for this step. This is the easiest done while laying the uprights on the floor. The rungs are to be placed into position as the illustration shows. The rungs should not extend out past the outside of the uprights once in place.

Once the rungs are in place, upright 'B' should then be placed on top of upright 'A' to hold the rungs in place. Use the provided hardware to fasten the two upright together. Tighten each of the fasteners being sure each rung is flush with the outside of the upright and are fitted properly as the rungs are oval in shape. After this step is complete, use the last rung and set it into place. Using two of the provided bolts and washers, tighten part 'C' into place.

Next, install the backboard to the wall. The hardware for this step is not included since we do not know your mounting surface. However, it is strongly advised that a professional be consulted for the installation of the backboard.

To properly mount the backboard, simply raise the upright and rung assembly to the wall and line up the backboard with the mounting holes being sure the bottom of the assembly is resting on the floor.

After the backboard is in place, the assembly can be raised and fastened to it. This step requires the use of 2-1/4-20 x 4" hex bolts and 2-flat washers.

It is recommended that the bolts on the whole assembly be checked periodically to be sure they are tight.

The Stall Bars are only meant for vertical weight support,

The Stall Bars measure 36" wide by 96" high.

