## Rolyan® Adjustable Outrigger Kit for MCP Flexion:

A832-1

Figures referred to in the instructions are on page 6.

### INDICATIONS

- To increase flexion of finger joints

### CONTRAINDICATIONS

- Not for use with fixed bony contractures
- Not for immediate treatment following fractures

### PRECAUTIONS

- **This splint is to be fitted initially by a healthcare professional who is familiar with the purpose for which it is prescribed. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.**
- **If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**
- **Be careful that the straps are not so tight that they interfere with circulation.**
- **When cutting metal end, hold or cover piece to be removed to prevent it from flying.**

### CARE AND CLEANING

All components can be cleaned with soap and lukewarm water. The straps can be washed with soap and water, but the adhesive portion should not be moistened. The splint and straps should be allowed to dry thoroughly before reapplication of the splint.

### INSTRUCTIONS FOR USE

1. Fabricate appropriate splint base.
2. Position outrigger on splint with extended portion to radial side and mark position (A).
3. Cut piece of low-temperature thermoplastic large enough to cover end sections of outrigger (B).
4. If splint base and overlapping piece are coated, remove coating by scraping surface or brushing it with solvent.
5. Heat end sections of outrigger with heat gun and embed into splint base, making sure ends do not poke through splinting material.
6. Heat overlapping piece of splinting material, make sure it is hot and dry, and place it over end sections of outrigger to bond them to splint base.
7. Attach straps to splint as appropriate.
8. With self-adhesive loop on surface of splint that will be toward skin, wrap D-ring rubber band attachment strap around proximal end of splint. Position strap so D-ring is on outside surface of splint (C).
9. Attach each rod adjuster loosely to outrigger by slipping slot in rod adjuster over outrigger wire and tightening set screws with Allen wrench (D).
10. Insert an extender rod through each rod adjuster. Adjust position of rod adjuster so extender rod lies parallel to joint being splinted. Use Allen wrench to tighten screw on rod adjuster to hold extender rod in place (E).
11. To determine final position of extender rod, use Allen wrench to provide a 90° angle relative to the finger (F).
12. Tighten the two set screws that hold rod adjuster on outrigger alternately, applying even tension to each screw.
13. Remove splint from patient and secure all set screws.
14. Using heavy-duty wire snips, cut off proximal end of extender rod. Extender rod should extend at least ½" (1.3 cm) beyond rod adjuster. Cover end with black tip protector.
15. Attach appropriate rubber bands or springs to proximal plastic D-ring (G).
16. With splint on extremity, position finger loop on finger and put monofilament into curled end of extender rod.
17. Attach spring or rubber band that is on proximal plastic D-ring through large hole in aluminum connecting piece on finger loop.
18. Adjust length of monofilament line to provide appropriate force of traction by sliding monofilament through two small holes on aluminum connecting piece (H). To increase amount of traction force, reduce length of monofilament. To reduce amount of traction force, increase length of monofilament.
19. A completed splint is shown in H.
20. For thumb flexion, attach a rod adjuster to curved ulnar portion of outrigger. Slide extender rod into rod adjuster and bend it to position finger loop to provide appropriate angle of pull (I).

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