Rev. 0 - 05/11\_107284

**Bauerfeind AG**  
Weissendorfer Strasse 5  
07937 Zeulenroda-Triebes  
Germany  
Phone +49 (0) 36628-66-35 00  
Fax +49 (0) 36628-66-39 99  
E-mail [info@bauerfeind.com](mailto:info@bauerfeind.com)

## GloboPed®

Vorfußentlastungsorthese „Postop“  
“Postop” forefoot-relieving orthosis  
L'orthèse de décharge de l'avant-pied « Postop »  
Ortesis de descarga del antepié „Postop“

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**GB** Dear Customer,

Please read through the instructions for use carefully and direct any questions to your physician, medical retailer, or our technical service department at the telephone number +44 (0) 121 446 53 53 (UK)/+1 800 423 34 05 (USA).

**Certification Information**

Bauerfeind AG is certified in accordance with DIN EN ISO 13485:2003 as well as DIN EN ISO 9001.

**Classification**

The **“Postop” forefoot-relieving orthosis** is a medical product as defined by the German Medical Product Law, classified in Class I (in accordance with Directive 93/92/EEC, Annex IX/III, 1.1).

**Service Number**

You can reach our Customer Service Department at the telephone number +44 (0) 121 446 53 53 (UK)/+1 800 423 34 05 (USA).

**Purpose**

The **“Postop” forefoot-relieving orthosis** is a high-quality orthopedic device for targeted relief for the area around the forefoot and front metatarsal area.

**Indications/Use**

- For early post-operation and/or post-traumatic recovery of the forefoot and front metatarsal area
- Treatment for diabetic foot ulcers in the forefoot area
- For rheumatic feet and algodystrophic syndromes

**Risks/Side Effects/Contraindications**

Mobility and balance may be compromised based on the purpose and desired function of the **forefoot-relieving orthosis**.

There are no known side effects on the entire body. The strap shoe should only be worn after consulting trained technical personnel.

There are no known hypersensitivities that cause illnesses. This device should only be worn after consulting your physician in the following clinical situations:

- Skin diseases/injuries in the body part being treated, particularly if inflammation appears. The same applies for accumulated scars with inflammation, redness, or overheating
- Acroparesthesia in the body part being treated
- Lymphoid circulation problems, unclear soft tissue inflammation around the device’s area of application

**Guaranty**

A guaranty is provided as stipulated by legal regulations. First-time application and fitting of the **forefoot-relieving orthosis** may be carried out by trained technical personnel only, otherwise no product liability will be assumed.

**Product Description**

The **“Postop” forefoot-relieving orthosis** is a prescribable product that may only be worn under a physician’s direction and only in accordance with these instructions of use and for the applications named.

The **“Postop” forefoot-relieving orthosis** may only be worn with a correctly installed protective frame that secures the forefoot against mechanical influence.

**Use**

This product is intended for the treatment of one individual patient and may not be re-used or used by others. No product liability will be assumed in the event of inappropriate use.

The **forefoot-relieving orthosis** should be put on while sitting and can be worn both on the left and right foot. Put on with the Velcro snugly on the foot, although not too tight in order to avoid pressure points on the skin respectively damages of arteries, and nerves.

Always wear a non-skid shoe with the same sole height on your healthy foot. If necessary, a physician may prescribe that the height of the sole / heel be changed by an orthopedic technician or orthopedic shoe technician.

**Safety Notice**

Any open wounds on the foot must be protected by a bandage or similar device. Do not ride a bicycle or drive a car or any other vehicle while wearing the orthosis.

Use crutches when wearing a **forefoot-relieving orthosis**.

Wearing the **“Postop” forefoot-relieving orthosis** is only recommended in hospitals and clinics or within your own home and not in public areas.

Please see your physician immediately should you notice any irregular changes in your body (i.e. increase in pain).

Do not alter the product in any way as it may affect the product’s functionality such that no product liability will be assumed.

**Technical Information**

Model	Size	Shoe Size	Article Number
Postop	S	33-34	3 662700 001
	M	35-38	3 662700 002
	L	39-42	3 662700 003
	XL	43-46	3 662700 004

**Accessories**

Model	Size	Article Number
Protective Frame, Black Plastic	S+M	3 662490 002
	L+XL	3 662490 003
Protective Frame, Aluminum	S-XL	3 662490 001

**Maintenance and Care**

Be sure to have your **“Postop” forefoot-relieving orthosis** checked regularly for functionality and closing functionality by an orthopedic technician or orthopedic shoe technician and exchanged if needed.

We recommend washing the orthosis separately at 40° C and by hand using a delicate detergent. We recommend closing the Velcro before washing in order to preserve its functionality and avoid damage. The orthosis should not be exposed to direct heat when drying (i.e. heating vents, sun light, etc.).

We have tested this product within the scope of our integrated quality management system. Should you nevertheless have a complaint, please contact our local technical consultant or our customer service department.

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Bauerfeind Middle East FZ-LLC
Bauerfeind USA, Inc.