## Clinically Proven Safe and Effective

Excerpts from a study conducted specifically and exclusively on Mother-To-Be® maternity support at the University of Oklahoma College of Medicine:

“No significant changes were encountered in the FHR baseline or beat-to-beat variability during placement of the support.”

“The maternal systolic, diastolic and mean arterial blood pressures were unaffected by the support. The right-sided and left-sided cardiac out-puts were unchanged during the monitoring periods.”

“This elastic/Velcro lumbar and abdominal support, available to relieve low backache, did not acutely affect the hemodynamics of the fetus and mother.”

“Each woman claimed that she felt more comfortable after being fitted with the support on the day of testing.”

“Each woman, when questioned two weeks later, reported improvement in back discomfort while sitting and standing.”

Reference: Journal of Reproductive Medicine: 1999; 44:1007-1011 “Low Backache During Pregnancy, Acute Hemodynamics Effects of a Lumbar Support” C.M. Beaty, M.D., V.J. Bhaktaram, M.D., E.F. Rayburn, M.D., et al

Manufactured by


P.O. Box 508 • 512 M Street  
Belleville, KS 66935  
(785) 527-5627

LRev 9/10 PL#4.504

CE

## To Fit Your *mother-to-be*® Maternity Support

Manufactured by

## To Wear • Façon de le porter • Instrucciones de uso

### 1. Putting On The Abdominal Support Pad

Place the Abdominal Support Pad (ASP) in a lifting position under the abdomen. Fasten it in place by wrapping its strap around your body and touching the Velcro hook to the ASP.

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### 3. Attaching The Belt

Grip the belt securely at both ends as shown. Stretch it forward and touch it to the ASP. You can adjust the tension and support you feel according to the amount you stretch the belt. To prevent skin irritation, be sure the Velcro hook on the ends of the belt is entirely on the ASP and doesn't touch your skin.

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### 2. Position The ASP

Place the ASP as low under your abdomen as possible but not so low as to interfere with your legs when you walk or sit. Slide the ASP up slightly until it doesn't interfere with your legs when lifted as shown. The ASP is now positioned correctly.

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
### 4. Fasten The Upper Tummy Strap

For an extra sense of security, lightly stretch the tummy strap over and touch it to the Velcro loop on the opposite side. The tummy strap is not essential but many prefer to have it. You can also use the strap to help hold the back support up while you are positioning it. The choice to use it or remove it is yours.

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### 5. Final Step (Optional)

Lay on your back with knees drawn up. Release the ends of the belt, stretch it again, re-attach it to the ASP, and stand up. This will lift the abdomen higher to achieve even greater relief from sagging weight and leg pain. Repeat this step each time you put on your support.

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>>> These may be purchased with or without heat moldable insert. May order insert for additional back support.

>>> French translation here. French translation here. French translation here. French translation here.

>>> Spanish translation here. Spanish translation here. Spanish translation here. Spanish translation here.

## Washing instructions

Hand wash in cold or warm water with mild detergent. Hang dry. Take care to cover the hook & loop so it doesn't tear the fabric as you handle it. Do not machine wash or dry.

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### Important Note:

Always consult your physician whenever you have questions about your health or your baby's.

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