





Polycentric Hinged Ulnar Deviation Splint

To reduce ulnar drift at the MP joints caused by rheumatoid arthritis, neurological impairment or trauma. Use following MP arthroplasty to maintain alignment after removal of post-operative splint. Recommended for ulnar drift that can be passively corrected with mild to moderate force.



The malleable finger stays may be opened to put the splint on or to take it off.



Gently contour each finger stay to create a half circle around the base of the finger. Avoid creating sharp curves in the finger stay. To prevent pressure in web spaces, stays may be rotated at their rivet attachments.



Apply wrist strap.

Splint shown on right hand.



Stevensville, MD 21666 U.S.A.

Made in the U.S.A.

ADJUSTMENTS

Edges of the Kydex shell may be flared with a heat gun. Avoid heating hinges or finger sections. Hold heat gun 1" to 2" away until plastic is pliable. Use glove or stockinette to handle hot Kydex. Allow to cool completely before fitting splint to hand.

TO ORDER

Measure the width across the 2nd to the 5th MP joints.

	Width at MP's	Right	Left
X-small	2 1/4" - 2 3/4"	P2002-R1	P2002-L1
Small	2 3/4" - 3"	P2002-R2	P2002-L2
Medium	3" - 3 3/8"	P2002-R3	P2002-L3
Large	3 3/8" - 3 5/8"	P2002-R4	P2002-L4

Order Toll Free **888.378.7763** (U.S only)

Local/Int. 410.604.6393

24 hour Fax 410.604.6398

Monday - Friday 8:30 - 5:30 pm EST

CARE INSTRUCTIONS



Wash with mild soap and water. Wrist strap can be hand washed and air dried.

Material: Kydex, Neoprene, Nylon Strap


©2005, 3-Point Products, Inc.

Authorized European Representative

CE
Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands

www.3pointproducts.com

With more severe ulnar deviation, choose the Radial Hinged Ulnar Deviation Splint.



Radial Hinged Ulnar Deviation Splint

See our full line of 3-Point splints

For night time rest and protection, choose the Comforter splint.



Comforter Splint