Instructions for use


Form Fit® Ankle Brace


Life Without Limitations®

.....	3
(EN) Instructions for use	5
(DE) Gebrauchsanweisung	6
(FR) Notice d'utilisation	7
(ES) Instrucciones para el uso	8
(IT) Istruzioni per l'uso	9
(NO) Bruksanvisning	10
(DA) Brugsanvisning	11
(SV) Bruksanvisning	12
(FI) Käyttöohjeet	13
(NL) Gebruiksaanwijzing	14
(PT) Instruções de Utilização	15
(ZH) 中文说明书	16

ENGLISH

USES

Mild to moderate sprains and strains; injury prevention, chronically weak ankles.

CAUTION

1. Patients should be instructed to inspect their foot, ankle, and this product before and after each use.
2. This product is designed to be worn over an athletic sock. If any problems occur with this product please contact your certified medical professional. Measure ankle circumference over an athletic sock around the base of the heel and across the ankle bone to compare with the sizing chart to ensure proper product fit.

PHOTO CAPTIONS

- A. Loop Material*
- B. Lace Hooks
- C. Lace Eyelets
- D. Tongue
- E. Figure-8 Straps*
- F. Finger Pull Tabs*
- G. Heel Opening

* Figure-8 version only

FITTING INSTRUCTIONS

1. Disconnect Figure-8 Straps (E)* and loosen laces
2. Insert foot into the brace toes first. Place heel securely into the Heel Opening (G). Adjust the Tongue (D) if necessary – making sure it is on top of the foot. (Fig.1)
3. Lace up the brace using either the Lace Hooks (B) (Fig.2) or the Lace Eyelets (C) (Fig.3).
4. For version with Figure-8 Straps: Wrap the inside (medial) Figure-8 Strap (E)* over the top of the laces and underneath the foot. Then fasten the end of the strap to the Loop Material (A) located on the inner side of the brace. Repeat for outside (lateral) Figure-8 Strap (E)* (Fig. 4).
5. Stand and apply full weight to the brace. Adjust Figure-8 Straps (E)* and laces to ensure a snug fit and maximum support (Fig. 5).

NOTE: During athletic activity the brace can be quickly adjusted to optimize comfort and support: use the Finger Pull Tabs (F) to adjust the tension on the Figure-8 Straps (E)*.

WASHING INSTRUCTIONS

1. This product may be cleaned using a mild detergent and 30°C (86°F) water. Remove shoe laces and hand wash separately. Fasten straps as if in use. Do not use bleach or strong detergents.
2. Lay brace out flat on a towel and air dry. Do not squeeze or wring out the brace. Do not place in dryer or near heater. Inspect brace carefully after each washing.

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Össur Iberia S.A.