Rolyan® AquaForm® Zippered Wrist Splint:

A175-001 through A175-003 and A175-010 through A175-012 and 0815-39-683, 691, 709, 659, 667, 675

Rolyan® AquaForm® Zippered Wrist & Thumb Spica Splint:

A175-100, 102, 104 and A175-110, 112, 114 and 0815-39-592, 600, 618, 626, 634, 642

Rolyan® AquaForm® Universal Zippered Wrist & Thumb Spica Splint:

5604-22 through 5604-25

Figures referred to in the instructions are on page 7.

Be sure to note which splint you are using and follow the instructions that apply to that splint. Unless otherwise stated, these instructions apply to both of the splints.

INDICATIONS

- May be used in the presence of arthritis, tendonitis and cumulative trauma injuries, and following sprains and strains

Rolyan AquaForm Zippered Wrist Splint

- For immobilizing the wrist
- May be used following tendon repair and to stabilize wrist fractures

Rolyan AquaForm Zippered Wrist and Thumb Spica Splint

- For immobilizing the IP, MCP, and CMC joints of the thumb in addition to the wrist (can be trimmed to allow IP movement)
- May be used for deQuervain's syndrome, dislocations, tenosynovitis of the CMC joint and gamekeeper's thumb, following healing of scaphoid fracture, and following MCP reconstruction

CONTRAINDICATIONS

- Not in the presence of severe edema and reflex sympathetic dystrophy
- May not be recommended for the acute phase of wrist fractures

PRECAUTIONS

- **These splints are to be fitted initially by a healthcare professional who is familiar with the purpose for which they are prescribed. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.**
- **If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**
- **Be sure the splint does not impinge the long extensor tendons of the thumb, the superficial radial nerve, and the dorsal branch of the ulnar nerve.**
- **After heating splinting material, always check its temperature before applying it to the patient.**

CARE AND CLEANING

- The splint will lose its shape in temperatures over 135°F (57°C). It should be kept away from sources of heat such as ovens, hot water, open flames, and sunny car windows.
- The splint may be cleaned with soap and lukewarm water. It should be allowed to air dry thoroughly before reapplication.

INSTRUCTIONS FOR USE

The zipper pull will be on the outside of the splint along the ulnar border when fabrication is complete.

Customizing the fit (to be done prior to heating splint)

1. Apply a stockinette gauntlet to cover the hand, wrist, and forearm. If softening the splint in a convection oven, use two layers of stockinette, with more over newly healed scars or skin grafts.
2. Pieces of 1/8" (3.2mm) low-tack **Polycushion Padding®** (Product No. A291-3) may be used to pad the radial and ulnar prominences. Cut the padding to fit and secure the pieces on the prominence(s) to be padded. They will be transferred to the splint after fabrication.
3. Position the splint on the patient to determine if trimming is necessary. If so, mark trim lines.
4. Trim the proximal edge as follows. (If the splint needs to be shortened by a small amount, the excess may be folded over [see instruction 14].)
 - Position the zipper pull at the distal border of the splint.
 - Heat the marked area in a heat pan of water (see instruction 5).
 - Cut along trim lines. Do not cut the zipper or fabric. Trim along the zipper about 1/4" (6.4 mm) from the fabric (A).
 - Lightly heat the splinting material that is sewn to the zipper and extends past the trim line. Also heat an equal distance distally along the zipper. Dry the area.
 - Apply solvent to the outside surface of the splint in these areas.
 - Fold the fabric and splinting material to the outside of the splint along the zipper and adhere the surfaces (B). Be sure to allow room for the zipper pull to move between the folded-over portions.

Molding the splint

Keep the stockinette liner smooth and wrinkle-free throughout splint fabrication and wearing.

5. Heat the 3/32" (2.4mm) **Aquaplast®-T** splint as follows:

Heat pan: Heat in water at a temperature of 160° to 170°F (70°–75°C) for one minute or until transparent. Remove from heat pan at an angle to let water drain off without stretching the brace and place on a flat surface. (The metal parts of the zipper will be hot.) Pat dry with a towel. The fabric around the zipper will be damp, but it will dry without special attention.

Convection oven: Heat without overlapping or folding on a baking sheet with a non-stick surface in an oven preheated to 195° to 205°F (90–95°C) until transparent (softening may take up to 10 minutes). Remove from baking sheet and place on a flat surface.

6. Position the wrist and forearm with the radial border up and the thumb in a functional position.
7. Be sure that the zipper pull is on the outside of the splint.

— **Wrist splint:** Place the thumb through the thumb hole (C).

— **Wrist and thumb spica splint:** Line up the longest extension on the distal edge of the splint with the thumb (D). Be sure that the splint edge follows the distal palmar crease.

For both splints, the body of the splint should drape on either side of the radial border of the forearm and the zipper sides should lie along the ulnar border.

8. **Wrist and thumb spica splint:** The distal portion of the splint extends on both sides of the thumb. Beginning with either side, wrap these two flaps around the thumb and overlap them (E).
9. For both splints, starting near the thumb, gently stretch the splint from the radial border around the forearm and hand until the zipper edges are touching. An even stretch of the material will result in a generally even pull on the perforations (F). The zipper should be parallel to the ulnar border. For comfort, the zipper should not be positioned directly over the ulnar prominence.
10. Close the zipper gently to ensure proper positioning of the extremity, the splint liner, and the splint (G).
11. Allow the splint to cool thoroughly and harden while on the extremity.

Modifying the fit

12. **Wrist and thumb spica splint:** The thumb portion may be trimmed to allow the IP joint free motion. If trimming is indicated, lightly heat only the area to be trimmed to avoid the possibility that the splint may lose its shape.
13. **Wrist and thumb spica splint:** The thumb piece may be held closed with straps to allow for fluctuating edema. It may also be bonded closed by following these steps:
 - Lightly heat the two edges where they overlap. Be sure the areas are dry.
 - Apply solvent to the areas.
 - Place the pieces together.
14. If desired, heat the proximal edge and flare it back. This technique can also be used to shorten the splint a small amount. For a stronger bond, apply solvent to the area. Gradually reduce the folded amount when nearing the zipper (H).

Wrist splint: The thumb hole area may also be flared.

Wrist and thumb spica splint: The distal edge of the thumb piece may be flared.

Completing the splint

15. Secure the optional padding pieces to the splint.
16. A completed wrist and thumb spica splint is shown in I.

(U.S. Patent 6,093,161)

Design Copyright Patterson Medical Company

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Patterson Medical, Ltd.**