


# LED Fiber Optic Curtain | 0815-55-614

- Includes LED fiber optic curtain, enclosure and a low-voltage LED light source
- Can be used for relaxation or de-escalation for individuals with dementia
- Great for theming and storytelling
- Ideal for developing cause and effect, color recognition skills
- Easily accessible for wheelchairs, standers and seating systems users thanks to the swiveling wall-mounted bracket
- Maintenance-free LED technology
- Operates on mains voltage transformed to low voltage



## **DIMENSIONS:**

Enclosure: 10" H x 28 ¾" W x 7" D  
Fiber Optic: 8' x 150 strands  
Controller: 13" diameter x 3 ¾" deep  
All dimensions are approximate.



1000 Remington Blvd., Suite 210, Bollingbrook, IL, 60440  
T: 1-800-323-5547 | W: [www.pattersonmedical.com](http://www.pattersonmedical.com)