

Interactive LED Bubble Tube | 0815-55-606

- Includes 40" LED Bubble Tube and 5-way wireless controller with 5 switches
- The user can control 4 colors (green, yellow, red and blue) visible from the bubble tube with the wireless controller
- The bubbles can be switched on and off with the wireless controller
- Ideal for distraction therapy and social interaction
- Great for developing color recognition, cause and effect, hand-to-eye coordination and gross motor skills
- Aids to the development and maintenance of dexterity skills through the use of the 5-way wireless controller
- The wireless controller is powered by integral batteries (supplied) and operates from a distance of up to 90ft
- Calming mode automatically starts if no button is activated for a period of 5 minutes
- Maintenance-free LED technology
- Operates on mains voltage transformed to low voltage



DIMENSIONS:

Tube Base: 16" diameter x 8" deep

Tube diameter: 6"

Controller: 13" diameter x 3 ¾" deep

All dimensions are approximate.



1000 Remington Blvd., Suite 210, Bollingbrook, IL, 60440

T: 1-800-323-5547 | W: www.pattersonmedical.com