

PROHEALTH

Dial Scales



MULTIPLE CAPACITIES AVAILABLE

Models D350 / D1130

*Well-designed and quick-and-easy-to-read, **Detecto's Pro-Health Series Scales** are well-suited for locker rooms, hotels, fitness clubs, pharmacies, nutritional counseling centers, offices, and private households. They feature fast weighing with the weight display immediately registering on the large 7-inch dial after stepping on the scale's slip-resistant platform.*



PRO HEALTH SCALES

- Both Pounds/Kilograms and Kilograms-Only Versions are Available.
- Easy-to-Read, Large 7-Inch Dial Registers Weight Quickly.
- Features Adjusting Scroll to Set the Zeropoint and Balance.
- Spacious Platform Accommodates Large Shoe Sizes.
- Low-Profile, Compact Design Requires Little Floor Space.

DETECTO MODEL D350 / D350K PROHEALTH SCALES

The precision movement of this sturdy low-profile scale with circular display guarantees exact results and the design with a lot of white and chrome gives the scale a fresh, timeless look. The powder-coated steel casing is built to last a lifetime. The spacious tread area on the

weighing platform is covered in easy-to-clean, black synthetic material for fast wipe-downs. The modern, elegant scale design will make a positive impression in any room. Durable features include reliable corrosion protection by galvanizing the frame and parts.



model D350

D350 SPECIFICATIONS	
CAPACITY	350 lb x 1 lb / 160 kg x 1 kg
DIMENSIONS (OVERALL):	13.25" W x 20.25" D x 4.5" H 336 mm W x 514 mm D x 114 mm H
DIMENSIONS (PLATFORM):	10.6" W x 10.6" D 269 mm W x 269 mm D
DIAL DIAMETER:	7" / 178 mm
SHIPPING WEIGHT:	17 lb / 8 kg



model D350K

D350K SPECIFICATIONS	
CAPACITY	160 kg x 1 kg
DIMENSIONS (OVERALL):	13.25" W x 20.25" D x 4.5" H 336 mm W x 514 mm D x 114 mm H
DIMENSIONS (PLATFORM):	10.6" W x 10.6" D 269 mm W x 269 mm D
DIAL DIAMETER:	7" / 178 mm
SHIPPING WEIGHT:	17 lb / 8 kg

DETECTO MODEL D1130 / D1130K PROHEALTH SCALES

A slim, low-profile design, economical pricing, and non-skid vinyl mat are just a few features of the Detecto D1130 and D1130K scales - proven performers for budget-minded customers. They feature a white, baked-enamel finish and personal, colorful indicators to mark one

person's or a whole family's weights. The large 7-inch-diameter dial allows for easy readability even for tall people. Designed for personal use, fitness centers, and health care facilities. Both pounds/kilograms and kilograms-only versions are available.



model D1130

D1130 SPECIFICATIONS	
CAPACITY	300 lb x 1 lb
DIMENSIONS (OVERALL):	12" W x 16.75" D x 2.5" H 305 mm W x 425 mm D x 64 mm H
DIMENSIONS (PLATFORM):	11.25" W x 11.75" D 286 mm W x 298 mm D
DIAL DIAMETER:	7" / 178 mm
SHIPPING WEIGHT:	9 lb / 4 kg



model D1130K

D1130K SPECIFICATIONS	
CAPACITY	130 kg x 500 g
DIMENSIONS (OVERALL):	12" W x 16.75" D x 2.5" H 305 mm W x 425 mm D x 64 mm H
DIMENSIONS (PLATFORM):	11.25" W x 11.75" D 286 mm W x 298 mm D
DIAL DIAMETER:	7" / 178 mm
SHIPPING WEIGHT:	9 lb / 4 kg

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

DETECTO
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Sold By: