

Interactive LED Fanlite | 0815-55-655

- Includes LED fanlite and microphone
- Eight modes available: three automatic and five voice-activated
- Features include sensitivity and speed controls
- Internal microphone and socket for external microphone
- Great for developing vocalization and basic life skills
- The unit is designed to be wall-mounted to ensure stability
- Maintenance-free LED technology
- Operates on mains voltage transformed to low voltage



DIMENSIONS:

40" W x 25" H x 5" D
Weight: 30 lbs



1000 Remington Blvd., Suite 210, Bollingbrook, IL, 60440
T: 1-800-323-5547 | W: www.pattersonmedical.com