

Sensory Corner | 0815-55-820

- Great starter kit
- Includes all 3 items that are always found in a sensory room: LED bubble tube, LED fiber optic and projector
- Great for developing basic life skills such as color recognition and switching skills
- Ideal to create a calming and relaxing environment
- All interactive products come complete with a wireless controller which allows the user to change the colors of the LED bubble tube and LED fiber optic
- The wireless controller is powered by integral batteries (supplied) and operates from a distance of up to 90ft
- Calming mode automatically starts if no button is activated for a period of 5 minutes
- Maintenance-free LED technology
- Operates on mains transformed to low voltage



Includes:

- 1 x Interactive LED bubble tube – 60" & 9-way wireless controller
- 1 x Curved bubble tube plinth - Medium 36" W x 36" L x 16" H
- 2 x Acrylic mirrors - 36" W x 48" H
- 1 x Bubble tube bracket
- 1 x Bubble tube pump and hose
- 1 x BCB - Bubble tube additive
- 1 x Interactive LED fiber optic bundle (includes Interactive LED light source, tails and 8-way controller)
- 1 x Projector bundle, including Solar 250 projector, wheel rotator, and 3 effect wheels (Deep, firework and a blank wheel)
- 1 x Spare lamp for projector Solar 250
- 2 x Calming sounds CDs
- 1 x Fitting instructions
- 1 x Manual



1000 Remington Blvd., Suite 210, Bollingbrook, IL, 60440
T: 1-800-323-5547 | W: www.pattersonmedical.com