



Classics

transport chair

- Comfortable 19" seating area and padded armrests
- Swing-away, removable foot rests
- Easily folds for storage and portability
- 7.5" wheels for easy maneuverability



A33677 0000



CAREX HEALTH BRANDS
SIOUX FALLS, SD 57101
MADE IN CHINA ©2012
www.carex.com

Join the Care Connection at
carex.com/cc for offers, product
information and care solutions.



Classics


transport chair



Comfortable 19" seating area

Weight capacity: 300 lbs.

How to attach footrests:



How to unfold backrest:



How to engage breaks:

Lift up to lock wheels and press down to unlock

