

CARPAL GEL WRIST SUPPORT

 Lohmann & Rauscher

Comfortable Wrist Support with Gel Pad for Carpal Tunnel Post-Op



- Soft, anatomically designed gel pad specially formulated with 100% medical grade mineral oil
 - Moisturizes skin on contact
 - Helps reduce post-op scarring from carpal tunnel surgery
 - Exudes mineral oil consistently throughout the pad's lifetime, thus providing continuous benefits
 - Ideal for keloid and hypertrophic scarring
- Dermatologist tested and non-allergenic, the gel pad will not support bacterial growth
- Gel pad helps cushion and protect sensitive scar tissue against shock, impact and vibration

- Removable, flexible plastic palmar stay provides semi-rigid support and additional stability when needed
- Breathable fabric with two-way stretch provides comfort and conformability to the hand/wrist area
- Contoured design and adjustable thumb strap allows for full use of fingers and thumb
- Easy and adjustable one-handed application and removal
- Latex free



Measure wrist circumference

S/M 5-7"
L/XL 7-9"

SIZE	ITEM NO.
S/M, Left	0816-04-271
S/M, Right	0816-04-297
L/XL, Left	0816-04-289
L/XL, Right	0816-04-305


PATTERSON
MEDICAL

www.pattersonmedical.com

Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada call 1-877-566-8611