

T655MS TREADMILL



The T655MS is rated for high workloads, long hours and is specifically designed for cardiac and physical rehabilitation. This treadmill comes loaded with unique features, programs and user amenities to help treat more patients more effectively.

- Bi-directional belt motion for retro walking and gait observations
- Decline to -3% grade for back and knee rehab
- Meets UL requirements for current leakage in a medical setting
- Full length medical side rails with pediatric lower sections for better patient stability
- 500 lb user rating, 24/7 operation with 5 year parts, 3 year labor warranty
- Eco-Glide auto lubrication system

USER WEIGHT	DIMENSIONS	RUNNING AREA	INCLINE RANGE	SPEED RANGE	DRIVE MOTOR
500 lbs / 227 kg	87.8 x 39.2 x 56.7 in / 220 x 95.5 x 144 cm	62 x 22 in / 158 x 56 cm	-3 - 15°	Reverse: 0.1 - 3 mph Forward: 0.1 - 12 mph	5.0 HP AC-Servo

Readouts

Calories, Speed, Time, Distance, Cal/HR, METs, Pace, Incline, Heart Rate, Weight Loss Zone, Cardio Zone, Active Zone

Workout Programs

Quick Start, Manual, Random (x 10,000), Interval (x3), Glute (x2) Hill (x3), 5K, 10K, Fat Burn, HRC (x3), Fitness Test (Bruce, Gerkin)

Features

- CSAFE Port
- Medical Handrails
- Optional Wired Controller
- Forward and Reverse Speed Function (-3 to 12 MPH)
- 0.1 SAFE-START™ Speed (Forward & Reverse)
- Contact and Telemetry HR
- Multi-ply 62 x 22 in. Belt
- Integrated 3-Speed Fan
- ECO-GLIDE™ Auto Lubrication System
- Optional SA Well +™
- Dual Headphone Jack with Built-In Volume Control
- Emergency Stop
- Centered Contact HR

Warranty

Full Commercial - no usage limit
Frame Lifetime
MyFlex Plus™ Lifetime
Drive Motor Lifetime

Parts 5 Years*
Wear Items 5 Years*
Labor 3 Years*
High Wear Items 90 days, no labor

Tech Specs

Power Requirements 15A / 120v / 60Hz - Dedicated Circuit
Plug Requirements NEMA / 5-15P
Unit Weight 390 lbs / 177 kg
Step Up Height 8.75 in / 22.23 cm
Micro Amp Leakage Components Medical Grade ≤ 300µA

Console



Tri-Color LED



Special Features



5.0 HP ECO-DRIVE™ Motor



Optional Wired Controller

*U.S. and Canada only.