



# HealthSmart® Digital Blood Pressure Monitors

## Clinical Accuracy Meets Style

### Standard Series Upper Arm

04-635-001 / 081698083

- + Memory: 2 user memory storage (60 readings total)
- + Display: date and time stamp, visual BP guide, average of last 3 readings and irregular heartbeat detection
- + Includes standard adult cuff (11.75"-16.50")



- ✓ 3-YEAR LIMITED WARRANTY
- ✓ CLINICAL ACCURACY

### Standard Series Universal Wrist

04-810-001 / 081698091

- + Memory: 2 user memory storage (120 readings total)
- + Display: date and time stamp, visual BP guide, average of last 3 readings and irregular heartbeat detection
- + Includes batteries, cuff (5.30"-8.50") and case



- ✓ 3-YEAR LIMITED WARRANTY
- ✓ CLINICAL ACCURACY

### Select Series Upper Arm

04-645-001 / 081698109

- + Memory: 2 user memory storage (120 readings total)
- + Display: backlit LCD screen, date and time stamp, visual BP guide, average of last 3 readings and irregular heartbeat detection
- + Upper arm unit includes standard adult cuff (11.75"-16.50"), AC adapter, batteries and zippered case



- ✓ 5-YEAR LIMITED WARRANTY
- ✓ CLINICAL ACCURACY

### Premium Series Universal Wrist

04-820-001 / 081698117

- + Memory: 2 user memory storage (120 readings total)
- + Display: high contrast back lit LCD screen, date and time stamp, visual BP guide, average of last 3 readings and irregular heartbeat detection
- + Audio readings in English or Spanish with volume control
- + Features heart spot technology, to indicate the optimum position of a user's arm for an accurate reading
- + Includes batteries, cuff (5.30"-8.50") and case



- ✓ LIFETIME LIMITED WARRANTY
- ✓ CLINICAL ACCURACY
- ✓ HEART SPOT TECHNOLOGY

### Premium Series Upper Arm

04-655-001 / 081698125

- + Memory: 2 user memory storage (120 readings total)
- + Display: high contrast back lit LCD screen, date and time stamp, visual BP guide, average of last 3 readings and irregular heartbeat detection
- + Audio readings in English or Spanish with volume control
- + Includes standard adult cuff (11.75"-16.50"), large adult cuff (16.50"-18.875"), AC adapter, batteries and zippered case



- ✓ LIFETIME LIMITED WARRANTY
- ✓ CLINICAL ACCURACY

### Accessories

04-240-001 – Std. Adult Cuff (11.75"– 16.50")

04-240-004 – Sm. Adult Cuff (8.625"– 14.125")

04-240-006 – Lg. Adult Cuff (16.50"– 18.875")

04-224-000 – AC Adapter UL

## Upper Arm Units

Specifications	Standard Series Upper Arm	Select Series Upper Arm	Premium Series Upper Arm
Item #	04-635-001	04-645-001	04-655-001
Retail Packaging	Yes	Yes	Yes
Packaging #	921-635-001	921-645-001	921-655-001
Packaging Dimensions	5" x 5.88" x 3.47"	5.12" x 9.19" x 5.03"	5.12" x 9.19" x 5.03"
Unit UPC	7-67056-63501-2	7-67056-64501-1	7-67056-65501-0
Unit Dimensions	3.9" x 2.9" x 5.28"	4.49" x 2.83" x 6.54"	4.49" x 2.83" x 6.54"
Unit Weight	14.30 oz.	14.99 oz.	17.21 oz.
Inner Carton Barcode	30767056635013	30767056645012	30767056655011
Inner Carton Quantity	3	3	3
Inner Carton Dimensions	11.20" x 6.30" x 5.20"	16.93" x 9.65" x 5.51"	15.38" x 9.65" x 5.51"
Inner Carton Weight	–	–	–
Master Carton Barcode	10767056635019	10767056645018	10767056655017
Master Carton Quantity	12	12	12
Master Carton Dimensions	12.20" x 11.41" x 14.37"	23.82" x 11.22" x 17.91"	23.82" x 11.22" x 17.91"
Master Carton Weight	–	–	–

## Wrist Units

Specifications	Standard Series Universal Wrist	Premium Series Universal Wrist
Item #	04-810-001	04-820-001
Retail Packaging	Yes	Yes
Packaging #	921-810-001	921-820-001
Packaging Dimensions	4" x 6.25" x 4.97"	4" x 6.25" x 4.97"
Unit UPC	7-67056-81001-3	7-67056-82001-2
Unit Dimensions	2.52" x 1.26" x 3.03"	2.52" x 1.26" x 3.03"
Unit Weight	3.88 oz.	4.13 oz.
Inner Carton Barcode	30767056810014	30767056820013
Inner Carton Quantity	4	4
Inner Carton Dimensions	10.55" x 6.69" x 8.58"	10.55" x 6.69" x 8.58"
Inner Carton Weight	–	–
Master Carton Barcode	10767056810010	10767056820019
Master Carton Quantity	24	24
Master Carton Dimensions	26.97" x 14.96" x 11.81"	26.97" x 14.96" x 11.81"
Master Carton Weight	–	–

## Accessories

Specifications	Std. Adult Cuff (11.75"– 16.50")	Sm. Adult Cuff (8.625"– 14.125")	Large Adult Cuff (16.50"– 18.875")	AC Adapter UL
Item #	04-240-001	04-240-004	04-240-006	04-224-000
Retail Packaging	Yes	Yes	Yes	Yes
Packaging #	921-240-001	921-240-004	921-240-006	921-240-000
Packaging Dimensions	–	–	–	–
Unit UPC	7-67056-24001-8	7-67056-24004-9	7-67056-24006-3	7-67056-22400-1
Unit Dimensions	–	–	–	–
Unit Weight	–	–	–	–
Inner Carton Barcode	30767056240019	30767056240040	30767056240064	7 30767056224002
Inner Carton Quantity	–	–	–	–
Inner Carton Dimensions	–	–	–	–
Inner Carton Weight	–	–	–	–
Master Carton Barcode	10767056240015	10767056240046	10767056240060	10767056224008
Master Carton Quantity	–	–	–	–
Master Carton Dimensions	–	–	–	–
Master Carton Weight	–	–	–	–

Specifications are subject to change without notice.