

Starfish™ Order Form

Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at **800-547-4333** or call in your order to our Pediatric Helpline at **800-499-0285**. Prices effective from January 1, 2014 through December 31, 2014.

Account No. _____
Date: _____
PO#: _____
Buyer: _____
Attention: _____

Ship to: _____
Name: _____
Address: _____
City: _____
Phone: () _____



Tumble Forms Starfish Bath Chair

Designed to reach new levels of functionality and ease of adjustability. Angle adjustable chair legs with articulating back and leg support, offer a range of height/angle positions to make bath time easier and more fun. Only one hand is needed to release the leg and seat locks located on both sides of the chair making bath time adjustments a breeze. The padded support bars give added positioning comfort.

The frame is aluminum, so it's easy to lift and carry. It folds flat for convenient storage or transport and includes nonslip foam pads for secure tub placement and protection. The seat features extra-soft, breathable mesh material that is removable for machine washing and comes in your choice of four fun colors. Available in four sizes to accommodate the youngest child through adolescence.



Back reclines to 0, 30, 60 and 90 degrees.

Adjustment locks can be released from either side with just one hand.

Leg support adjusts to -60, -30, 0, 30 and 60 degrees.

Folds flat making it easy to carry and store.



Set rear leg to -30° to tilt the seat for a level of versatility not found on other bath chairs.

Starfish Bath Chair

Size 1	Size 2	Size 3	Size 4	
___ 9268-76-02	___ 9268-77-02	___ 9268-78-02	___ 9268-79-02	Ocean Blue
___ 9268-76-04	___ 9268-77-04	___ 9268-78-04	___ 9268-79-04	Lobster Red
___ 9268-76-05	___ 9268-77-05	___ 9268-78-05	___ 9268-79-05	Puffer Fish Purple
___ 9268-76-06	___ 9268-77-06	___ 9268-78-06	___ 9268-79-06	Jellyfish Pink
\$479.70	\$492.30	\$523.85	\$536.50	



Lobster Red



Ocean Blue



Jellyfish Pink



Puffer Fish Purple

Accessories for Starfish Bath Chair

Abduction Strap

Size 1 & 2	Size 3 & 4
___ 9215-24-2102	___ 9215-27
___ 9215-24-2104	___ 9215-29
___ 9215-24-2105	___ 9215-30
___ 9215-24-2106	___ 9215-31
\$54.25	\$56.75

Head Laterals

One Size	
___ 9215-24-3102	Ocean Blue
___ 9215-24-3104	Lobster Red
___ 9215-24-3105	Puffer Fish Purple
___ 9215-24-3106	Jellyfish Pink
\$54.25	

Starfish Bath Chair Replacement Covers

Size 1	Size 2	Size 3	Size 4	
___ 9215-24-0102	___ 9215-24-0202	___ 9215-24-0302	___ 9215-24-0402	Ocean Blue
___ 9215-24-0104	___ 9215-24-0204	___ 9215-24-0304	___ 9215-24-0404	Lobster Red
___ 9215-24-0105	___ 9215-24-0205	___ 9215-24-0305	___ 9215-24-0405	Puffer Fish Purple
___ 9215-24-0106	___ 9215-24-0206	___ 9215-24-0306	___ 9215-24-0406	Jellyfish Pink
\$89.60	\$89.60	\$98.45	\$98.45	

Starfish Bath Chair User Weight & Height			
Size	Max Weight (kg)	Min. Height (mm)	Max Height (mm)
1	159 lbs. (72)	31.5" (800)	39" (1000)
2	159 lbs. (72)	35" (900)	47" (1200)
3	159 lbs. (72)	43" (1100)	55" (1400)
4	159 lbs. (72)	51" (1300)	65" (1650)

Starfish Bath Chair Specifications				
Size (mm)	1	2	3	4
Head Support	6" (150)	6" (150)	8.5" (215)	8.5" (215)
Back Support Height	14" (350)	16" (400)	18" (460)	22" (550)
Seat Depth	9" (225)	11" (275)	13" (335)	16" (400)
Seat Width	16" (400)	16" (400)	16" (400)	16" (400)
Leg Support	10" (260)	10" (260)	13" (340)	13" (340)
Overall Length	39" (975)	42" (1,050)	52.5" (1,312.5)	52.5" (1,312.5)
Seat Height	2-17" (50-425)	2-17" (50-425)	2-17" (50-425)	2-17" (50-425)



Shower Base for Starfish Bath Chair. From the tub to the shower. The optional shower base easily attaches to the Starfish bath chair to bring the child to the ideal height within a strong, stable, maneuverable base. 14" from floor to base. Width: 21". Length: 34"–48". Seat height from floor: 20"–28".
9215-25 \$555.45