

Portable and Innovative

The Portable Pro Professional Grade is a strong and durable clinical grade ultrasound in a portable package. This device has 3 intensity settings and comes with a carrying case, gel and an AC adapter.

3 Power Intensity Modes and 3 Timer Modes



566193

Technical Specifications

Working Frequency:	1.0MHz+-10%
Max Value Power Output:	9.60W+-20% (Duty Cycle 100%)
Power Output:	L 2.88W, M: 3.84W, H: 4.80W
Modulate Frequency:	50Hz+-10%
Modulate Duty Cycle:	Adjustable, L (30%), M (40%), H (50%)
Effective Radiation Area:	4.0cm ² +-20%
Maximum Effective Sound Intensity:	2.4W/cm ² +-20% (Duty Cycle: 100%)
Max RBN:	5.0
Beam Type:	Collimated
Material of US Head:	Aluminum
Working Voltage:	15VDC
Working Current:	Working current less than 1.0A, Static current less 80mA
Working Time:	Adjustable 5 min, 10 min, 15 min
Max Controlled Temperature of Ultrasound:	42+-2° C (36° F)
Size:	202mm(L)x49mm(W)x0mm(H) (8"x2"x2¾")
Safety Style:	Class IIa, BF-type



1000 Remington Blvd, Suite 210
Bolingbrook, IL 60440

www.pattersonmedical.com
1-800-323-5547

