

BARRINGTON DIAGNOSTICS

BV-11-776NL **Wall Aneroid**

- Large, easy-to-read, 6" faceplate with white luminescent numerals
- Adult size cuff, fits arm circumference range 23-40 cm (10 3/4" to 15 3/4")
- Square plastic casing encloses a no stop pin, 300 mmHg gauge
- Calibration tool allows for on-site calibration
- LIFETIME calibration warranty +/- 3mmHg
- 8' Latex Free coiled tubing (black)
- Deluxe air release valve
- Easy viewing with swivel feature



New look Latex-Free Tubing (end view)



BV-11-786NL **Floor Aneroid**

- Includes features of model BV-11-776NL plus...
- Adjustable floor stand (42.5" - 55")
- Easy assembly
- * Stand can be purchased separately

Barrington Diagnostics
28W206 Commercial Ave. Unit B
Barrington, IL 60010
Ph: 847-382-7770 • Fax: 847-382-7790
Toll Free: 888-822-8293