


Soft Transfer Swivel

Item#5528-93



DESCRIPTION:

The soft Transfer Swivel conforms to the shape of a seat and swivels easily. When back and/or hip pain is present, the 16" diameter cushion provides ease of getting in and out of an automobile. The top is a soft cotton material that is filled with foam for comfort. The bottom is covered with slip-proof material for added safety. It is lightweight and portable. Weight limit: 350 lbs.



CARE:

To wash, separate halves and remove metal disc in middle. Machine wash using gentle cycle or by hand and air dry.



Sammons Preston
APattersonMedicalCompany
1000RemingtonBlvd.,Suite210
Bolingbrook,IL60440
Phone: 800-323-5547
DirectPhone:630-378-6000
Fax:800-547-4333
DirectFax:630-378-6310



Sammons Preston
APattersonMedicalCompany
W68N158EvergreenBlvd.
Cedarburg,WI53012



Authorized European Representative
Homecraft Rolyan
APattersonMedicalCompany
Nunn Brook Road, Huthwaite
Sutton-In-Ashfield, Nottinghamshire NG17
2HU
United Kingdom
Phone: 08702423305
Fax: 08701660211

International
Phone: +441623759146
Fax: +441623726886

©November 2010
Patterson Medical Products, Inc.