

Impulse 3000 with Timer

Size	3.8" x 2.60" x 1" / (97 mm x 62 mm x 25 mm)
Weight	3.1 oz. (87 g)
Power Source	9V Battery, E-block, type 6F22
Channels	Dual
Waveform	Asymmetrical, biphasic square-wave
Amplitude	0 – 98 (mA)
Frequency	2 - 150 (Hz or pps)
Pulse Width	50 - 250 microseconds (μ s)
Output	Constant current
Continuous	Continuous stimulation
Burst	8 pulses per burst / 2 bursts per second
Modulation	Decrease/increase of set Pulse Width of 50% over a five-second period
Timer	15, 30, 60 minutes
Number of Electrodes	4
Number of Lead wires	2
Tolerances	+/- 10%
Data recorded across a 500 OHM load resistance	