



APS Dry Needles are used by practitioners worldwide. Developed in cooperation with international experts specifically for **Dry Needling**, Physical Therapists no longer have to improvise with acupuncture needles!

Safe, high quality APS Needles improve the function of Dry Needling. The result is better outcomes and greater patient satisfaction. Whereas acupuncture needles were designed to be superficially inserted only once, the **Dry Needling** technique requires multiple deep insertions and manipulations on each of the targeted trigger points.

APS Needle Characteristics

- ◆ Headless needle for better precision when puncturing.
- ◆ Longer, thicker, stiffer handle improves handling and control throughout the treatment.
- ◆ High strength, patented Surgical Stainless Steel formula prevents bending for greater accuracy and less pain during treatment.
- ◆ Special conical sharpening improves needle penetration for multiple insertions.
- ◆ Triple lubricant coating maintains lubrication during the entire procedure.
- ◆ Triple polishing removes surface irregularities achieving a smooth painless needle.

Benefits

- ◆ A needle specially designed for Dry Needling.
- ◆ Perfect grip that ensures total control and precision.
- ◆ Easy penetration and better handling of the needle from beginning to the end of treatment.
- ◆ Smooth insertion minimize tissue damage.
- ◆ Reduction of pain during treatment.

Features

- ◆ Easy open packaging.
- ◆ Thirteen Sizes available including 100mm length.
- ◆ APS Dry Needling Needle sizes are differentiated with color tips for easy identification.