



3M™ Tegaderm™ Absorbent Clear Acrylic Dressing

A unique transparent wound dressing with fluid handling capacity.

Versatile—ideal for skin tears, pressure ulcers, abrasions and donor sites.

- Transparent dressing allows for wound monitoring without changing the dressing
- Clear design takes the guesswork out of application over the wound
- Novel acrylic polymer pad technology designed to handle low to moderate wound drainage



Order #	3M #	Dressing Pad Size	Overall Dressing Size	Unit of Measure
0815-74-821	90800	1½ in. x 2¼ in.	3 in. x 3¾ in.	5 Dressings/Box
0815-74-839	90801	2¾ in. x 3 in.	4¾ in. x 5 in.	5 Dressings/Box
0815-74-847	90803	3¾ in. x 4¼ in.	5½ in. x 6¼ in.	5 Dressings/Box
0815-74-854	90807	4½ in. x 5 in.	6½ in. x 7½ in.	5 Dressings/Box



www.pattersonmedical.com

To order call:
1-800-323-5547
 In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
 In Canada fax 1-877-566-8611