



Thera-Band Roller Massager+ Fact Sheet

General Description

The Green Thera-Band® Roller Massager+ is an innovative tool for myofascial release and deep tissue massage. Its unique patent-pending ridged design supports both superficial and deep tissue mobilization while providing a massage-like experience. Meanwhile, the roller massager handles were designed to facilitate trigger point release. Use of the Roller Massager+ can help increase blood flow and circulation in targeted areas, while also helping to increase muscle flexibility and range of motion.

Available in a Standard version that's great for use in-clinic and at home, plus a Portable version with patent-pending retractable handles that's perfect for travel and on-the-go applications. All Thera-Band Roller Massagers are individually packaged with detailed instructions featuring general techniques and advanced exercises divided by muscle region. They are constructed with durable **latex-free** thermoplastic material for extended life and easy cleaning. Thera-Band Roller Massagers can be used in-clinic, as a part of a home exercise program or patient/client self-initiated wellness routine.

Features & Benefits

- **Patent-pending ridged design**
 - Supports both superficial and deep tissue mobilization with unique ridged profile
- **Portable version**
 - Features patent-pending retractable handles for easy travel
- **Highly durable thermoplastic material**
 - Offers extended life and easy cleaning
 - Latex-Free
- **Multi-use handles**
 - Designed to support trigger point release
- **Detailed instructions**
 - Enable users to get the most out of their purchase in-clinic and at home
 - Color poster illustrates general techniques and advanced exercises by muscle region
 - Additional videos available online at: info.Thera-BandAcademy.com/RollerMassager

Packaging Details

Part #	Thera-Band Roller Massager+
5661-48	Standard, Green (Ridged)
5661-49	Portable, Green (Ridged)

Care & Maintenance

Store your Thera-Band Roller Massager at room temperature in a cool, dry place, outside of direct sunlight.

To clean the roller massager, we recommend using a mild detergent and rinse thoroughly in clean tap water. Allow roller massager to air dry out of direct sunlight.



Standard



Portable



1000 Remington Blvd., Suite 210
Bolingbrook, IL 60440
1-800-323-5547
www.pattersonmedical.com