

## 3pp™ Step Up™ Splint



Restore PIP and/or DIP extension.



3pp™ Side Step™ Splint




3pp™ Step Down™ Splint


**3point**  
products



## 3pp™ Step Up™ Splint




1




2



3



4



5



6

Watch for our 3pp videos at [www.youtube.com/3pointproducts](http://www.youtube.com/3pointproducts).

### Sizing

Measure length from web space to fingertip and width of proximal phalanx.

#### Web Space to tip

#### Width at Prox Phalanx

Sm	1 1/2" - 2"	(3.8 - 5 cm)	1/2" - 5/8"	(1.3 - 1.6 cm)
Med	1 3/4" - 3 1/8"	(4.4 - 8.3 cm)	5/8" - 7/8"	(1.6 - 2.2 cm)
Lrg	3" - 3 1/2"	(7.6 - 8.9 cm)	7/8" - 1 1/8"	(2.2 - 2.9 cm)

### To Order

Sm P1200-2

Med P1200-3


Large P1200-4

Order Toll Free **888.378.7763** (U.S only)

Local/Int. 410.604.6393

24 hour Fax 410.604.6398

Mon. - Thurs. 8:30 - 5:30 pm, Fri. 8:30 - 5:00 pm EST



Hand or machine wash straps and pads. Air dry.

#### Material

Steel wire, nylon and Lycra®, polyethylene foam. Latex free.

Made in Mexico

©2009, 3-Point Products, Inc.