

3pp™ Side Step™ Splint



Reduce lateral deformities of the PIP and DIP joints.



3pp™ Step Up™ Splint



3pp™ Step Down™ Splint

3point
products




3pp™ Side Step™ Splint



Slide up for PIP



Slide down for DIP



Apply hook

1



2



3

Watch for our 3pp videos at www.youtube.com/3pointproducts.

Sizing

Measure length from web space to fingertip and width of proximal phalanx.

Web Space to tip

Sm 1 1/2" - 2" (3.8 - 5 cm)

Med 1 3/4" - 3 1/8" (4.4 - 8.3 cm)

Lrg 3" - 3 1/2" (7.6 - 8.9 cm)

Width at Prox Phalanx

1/2" - 5/8" (1.3 - 1.6 cm)

5/8" - 7/8" (1.6 - 2.2 cm)

7/8" - 1 1/8" (2.2 - 2.9 cm)

To Order

Sm P1201-2

Med P1201-3

Large P1201-4

Order Toll Free 888.378.7763 (U.S only)

Local/Int. 410.604.6393

24 hour Fax 410.604.6398

Mon. - Thurs. 8:30 - 5:30 pm, Fri. 8:30 - 5:00 pm EST



Hand or machine wash straps and pads. Air dry.

Material

Steel wire, nylon and Lycra®, polyethylene foam. Latex free.

Made in Mexico

©2009, 3-Point Products, Inc.