

7910 Elastic Knee



Small 081693712
 Medium 081693720
 Large 081693738
 XLarge 081693746

Anatomically designed with no seams so the elastic support will comfortably contour to the shape of your knee. It provides even compression and support for weak and injured knees, while minimizing bunching.

7920 Elastic Knee Stabilizer



Small 081693753
 Medium 081693761
 Large 081693779
 XLarge 081693787

Provides compression & support for weak and injured knees. Anatomically designed for comfort and fit. Side stabilizers provide increased medial and lateral support to weak & injured knees.

Measure slightly bent underneath knee cap	
Small	12 - 13 1/2"
Medium	13 1/2 - 15"
Large	15 - 16 1/2"
XLarge	16 1/2 - 18"

7900 Elastic Ankle



Small 081693928
 Medium 081693936
 Large 081693944
 XLarge 081693951

Anatomically designed with no seams so the elastic support will comfortably contour to the shape of your ankle. Provides even compression & support for weak and injured ankles, while minimizing bunching.

7905 Elastic Ankle Wrap



S/M 081693969
 L/XL 081693977

Fits snugly around your ankle to provide compression and support. An adjustable closure strap allows you to apply additional compression either medially or laterally when needed.

Measure around ankle joint	
Small	6 1/4 - 8 1/4"
Medium	8 1/4 - 10 1/4"
Large	10 1/4 - 12 1/4"
XLarge	12 1/4 - 14 1/4"

7940 Elastic Back Stabilizer



Small 081693985
 Medium 081693993
 Large 081694009
 XLarge 081694017

Provides compression & support for general lower back dysfunctions. Built-in flexible stabilizers provide additional support to the lower back. Two adjustable side straps for extra compression.

Measure circumference of the waist	
Small	25 1/2 - 29 1/2"
Medium	29 1/2 - 33 1/2"
Large	33 1/2 - 37 1/2"
XLarge	37 1/2 - 41 1/2"

Call toll free:
1-800-323-5547

Fax toll free:
1-800-547-4333



PATTERSON
 MEDICAL

www.pattersonmedical.com