

# CaligaLoc®

Stabilizing orthosis for partial immobilization of the ankle.



Support and postural splint



## Indications

- Conservative treatment of torn ankle ligaments
- Postoperative protection after ligamental suturing/reconstruction
- Temporary stabilization for post-traumatic tarsal sinus syndrome and decompensated instability of the subtalar joint
- Permanent stabilization for chronic instability of the upper and/or lower ankle where surgery is contraindicated

## Features

- Logical outcome of the development of the proven MHH ankle splint.
- The stability element is anatomically contoured and can be reshaped if necessary.
- Can be worn around the clock, as a support splint in normal shoes during the day and as a postural splint, without shoes, at night.

## Benefits

- Stabilizes the ankle by counteracting supination.
- Prevents abnormal load situations in the upper and, to a certain extent, lower ankle.
- The integrated pronating raised heel section takes the strain off the lateral malleolar ligaments.

U.S. Shoe Size	1	2	3	Order number
Men	2½ - 6	7 - 9	10½ - 14	12013003080+0°
Women	3½ - 7½	8½ - 10½	12 - 15½	

+6 = right, +7 = left, ° = size