




EXOS

SIZING FOR BRACES



MEASURE
CIRCUMFERENCE OF HAND
OR USE
EXOS SIZING TAPE
(#800-SIZING TAPE)
AS SHOWN

BRACE SIZE	MEASUREMENT
3XSML	5 $\frac{3}{4}$ " - < 6 $\frac{1}{2}$ "
2XSML	6 $\frac{1}{2}$ " - < 7 $\frac{1}{4}$ "
XSML	7 $\frac{1}{4}$ " - < 8"
SML	8" - < 9"
MED	9" - < 10"
LG	10" - < 11"
XLG	11"