




LOWER EXTREMITY SIZING FORM

To order, please use the ReadyWrap Lower Extremity Order Form



New

New

New

ReadyWrap Thigh		Length		Girths	
Item Code	Item Desc.	Thigh	Above Knee	Groin	
RW-LE-EG-S25	Average, Small	up to 25 cm	38 – 48 cm	48 – 63 cm	
RW-LE-EG-M25	Average, Medium	up to 25 cm	45 – 58 cm	58 – 73 cm	
RW-LE-EG-L25	Average, Large	up to 25 cm	52 – 68 cm	68 – 83 cm	
RW-LE-EG-XL25	Average, X-Large	up to 25 cm	59 – 78 cm	78 – 93 cm	
RW-LE-EG-XXL25	Average, XX-Large	up to 25 cm	66 – 88 cm	88 – 103 cm	
RW-LE-EG-S30	Tall, Small	up to 30 cm	38 – 48 cm	48 – 63 cm	
RW-LE-EG-M30	Tall, Medium	up to 30 cm	45 – 58 cm	58 – 73 cm	
RW-LE-EG-L30	Tall, Large	up to 30 cm	52 – 68 cm	68 – 83 cm	
RW-LE-EG-XL30	Tall, X-Large	up to 30 cm	59 – 78 cm	78 – 93 cm	
RW-LE-EG-XXL30	Tall, XX-Large	up to 30 cm	66 – 88 cm	88 – 103 cm	

ReadyWrap Knee		Girths	
Item Code	Item Desc.	Below Knee	Above Knee
RW-LE-DE-S	Small	28 – 38 cm	38 – 48 cm
RW-LE-DE-M	Medium	33 – 43 cm	45 – 58 cm
RW-LE-DE-L	Large	38 – 48 cm	52 – 68 cm
RW-LE-DE-XL	X-Large	43 – 58 cm	59 – 78 cm
RW-LE-DE-XXL	XX-Large	53 – 68 cm	66 – 88 cm

ReadyWrap Calf		Length		Girths	
Item Code	Item Desc.	Calf	Ankle	Calf	
RW-LE-BD-S30	Average, Small	up to 30 cm	18 – 23 cm	28 – 38 cm	
RW-LE-BD-M30	Average, Medium	up to 30 cm	20 – 28 cm	33 – 43 cm	
RW-LE-BD-L30	Average, Large	up to 30 cm	25 – 33 cm	38 – 48 cm	
RW-LE-BD-XL30	Average, X-Large	up to 30 cm	30 – 38 cm	43 – 58 cm	
RW-LE-BD-XXL30	Average, XX-Large	up to 30 cm	35 – 43 cm	53 – 68 cm	
RW-LE-BD-S35	Tall, Small	up to 35 cm	18 – 23 cm	28 – 38 cm	
RW-LE-BD-M35	Tall, Medium	up to 35 cm	20 – 28 cm	33 – 43 cm	
RW-LE-BD-L35	Tall, Large	up to 35 cm	25 – 33 cm	38 – 48 cm	
RW-LE-BD-XL35	Tall, X-Large	up to 35 cm	30 – 38 cm	43 – 58 cm	
RW-LE-BD-XXL35	Tall, XX-Large	up to 35 cm	35 – 43 cm	53 – 68 cm	

ReadyWrap Foot		Length		Girths	
Item Code	Item Desc.	Foot	Y	Ankle	
RW-LE-AB-L/R-S	Left / Right, Small	up to 21 cm	28 – 33 cm	18 – 23 cm	
RW-LE-AB-L/R-M	Left / Right, Medium	up to 23 cm	30 – 38 cm	20 – 28 cm	
RW-LE-AB-L/R-L	Left / Right, Large	up to 26 cm	35 – 43 cm	25 – 33 cm	
RW-LE-AB-L/R-XL	Left / Right, X-Large	up to 29 cm	40 – 48 cm	30 – 38 cm	

ReadyWrap Toe		Girth	
Item Code	Item Desc.	MTP	
RW-LE-AA-L/R-S	Left / Right, Small	18.0 – 20.4 cm	
RW-LE-AA-L/R-M	Left / Right, Medium	20.5 – 23.5 cm	
RW-LE-AA-L/R-L	Left / Right, Large	23.6 – 26.0 cm	
RW-LE-AA-L/R-XL	Left / Right, X-Large	26.1 – 28.5 cm	
RW-LE-AA-L/R-XXL	Left / Right, XX-Large	28.6 – 31.0 cm	

Terms & Conditions: Net 30 terms available for customers with approved credit applications. Other orders must be prepaid by check, Mastercard, Visa, or Discover.
Returns & Cancellations: Please inspect your order immediately upon receipt. Defective products may be returned or exchanged within two weeks of delivery. Products returned which are not defective are subject to round-trip shipping charges and a restocking fee of 15%. Prices and policies are subject to change without prior notification.