



Landing Strip® Bedside Fall Injury Prevention

Indications for Use:

Position the Landing Strip bedside. It is designed to remain on the floor **at all times**, especially if the patient is at risk for falls or has difficulty transferring. Just like an airbag in your car, the Landing Strip is designed to remain in place at all times. You never know when an accident will occur.

- Can be used with low or standard height frames.
- Compatible use for those who are ambulatory, those who use walkers, or those in wheelchairs.
- The lighter color is recommended for dementia patients. Some may perceive dark, rich colored mats as an obstacle or a hole in floor.

Suitable for people:

- Anyone who has experienced a previous fall
- Assess medical condition (e.g. dizziness, poor hydration or nutrition)
- Mental Status (confused or agitated)
- Disorientation (forgets limitations)
- Impaired or weak gait
- Difficult or assist required transferring
- IV therapy heparin lock
- Bladder control issues (may attempt unassisted ambulation)



Available in gray & burgundy.
Sizes available: 24" x 70" x 1" & 36" x 70" x 1".



Cleaning Instructions: Wipe clean with a disinfectant solution. Keep non-skid bottom clean and allow to dry before placing it bedside.



www.pattersonmedical.com

Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada Fax 1-877-566-8611