

# circaid<sup>®</sup> reduction kit

An innovative system for acute phase lymphedema compression.

## Indications:

- Early phase decongestion treatment of lymphedema.

## Features:

- Inelastic compression.
- Breathable fabric that improves comfort and compliance.
- Juxtaposed band system for a perfect fit.
- Easily customized on a per-patient basis.
- Easy handling when donning and doffing.
- Built-In Tension System supports a repeatable level of tension.
- Adjustable to your specific patient's needs.

## System Includes:

- Upper leg system includes: 1 reduction upper leg component, 1 reduction shelf strap, 6 hook and loop stays, 1 Built-In Tension system guide card, 1 paper measuring tape, and 1 direction for use.
- Knee system includes: 1 reduction knee component, 2 hook and loop knee spines, 1 paper measuring tape, 1 direction for use, and 1 reduction shelf strap.
- Lower leg system includes: 1 reduction lower leg component, 1 pair reduction kit undersleeve leg, 1 reduction shelf strap, 6 hook and loop stays, 1 Built-In Tension system guide card, 1 paper measuring tape, and 1 direction for use.
- Arm system includes: 1 reduction arm component, 1 pair undersleeve arm, 1 customizable hand wrap, 1 Built-In Tension system guide card, and 1 direction for use.
- Foot options are not included with circaid<sup>®</sup> reduction kits and must be ordered separately.



circaid<sup>®</sup> reduction kit whole leg with undersock



circaid<sup>®</sup> reduction kit arm with undersleeve

Now available from




[www.pattersonmedical.com](http://www.pattersonmedical.com)

Call toll free:  
**1-800-323-5547**

Fax toll free:  
**1-800-547-4333**

# Now available from PATTERSON MEDICAL

medi



<b>circaid reduction kit upper leg</b>			<b>Thigh Length</b>		
			Short <37cm	Standard 37-41cm	Long >41cm
<b>Circumference</b>	Regular	Single	0816-85-767	0816-85-866	0816-85-882
		4-pack	0816-85-775	0816-85-874	0816-85-908
	Wide	Single	0816-85-932	0816-85-957	0816-85-973
		4-pack	0816-85-940	0816-85-965	0816-85-981

<b>circaid reduction kit lower leg</b>			<b>Calf Length</b>		
			Short <32cm	Standard 32-36cm	Long >36cm
<b>Circumference</b>	Regular	Single	0816-85-494	0816-85-510	0816-85-536
		4-pack	0816-85-502	0816-85-528	0816-85-544
	Wide	Single	0816-85-551	0816-85-577	0816-85-593
		4-pack	0816-85-569	0816-85-585	0816-85-601


<b>circaid reduction kit knee</b>		<b>Knee Length</b>
		30cm
<b>Circumference</b>	Single	0816-85-734
	4-pack	0816-85-742

reduction knee kit is only available in one length circumference.

<b>Leg Measurement Guide</b>	
T	Waist
H	Widest part of pelvis
K	Gluteal fold
G	Groin
F	Mid-thigh
E2	15cm above E
E1	5cm above E
E	Center of patella
D1	Knee crease
D	Slightly below knee
C	Widest part of calf
B1	Between ankle and widest part of calf
B	Narrowest part of ankle above malleolus
A1	Middle of foot
A	Ball of foot
Y	Diagonally around heel over widest part of top of ankle
Z	Heel to base of great toe
mall	Center of malleolus
ℓ = length	
c = circumference	

<b>circaid reduction kit arm</b>			<b>Arm Length</b>	
			Standard <45cm	Long >45cm
<b>Circumference</b>	Regular	Single	0816-91-245	0816-91-260
		4-pack	0816-91-252	0816-91-278
	Wide	Single	0816-91-286	0816-86-005
		4-pack	0816-85-999	0816-86-013

<b>circaid reduction kit accessory items</b>	
reduction kit undersleeve leg pair	0816-85-692
reduction kit undersleeve leg pair 4-pack	0816-85-700
reduction kit undersleeve arm standard 8-pack	0816-86-039
reduction kit undersleeve arm wide 1-piece	0816-86-047
reduction kit undersleeve arm wide 8-pack	0816-86-054
circaid pac band xlarge	0816-85-650
circaid pac band xlarge 8-pack	0816-85-668
circaid customizable hand wrap	0816-86-021
circaid compression anklet pair standard 8-pack	0816-85-676
circaid compression anklet pair large 8-pack	0816-85-684
reduction kit extension short 30cm	0816-85-619
reduction kit extension standard 35cm	0816-85-627
reduction kit extension long 40cm	0816-85-635
circaid single band ankle foot wrap 8-pack	0816-85-643
reduction kit shelf strap	0816-85-718
reduction kit shelf strap 8-pack	0816-85-726
reduction kit knee spine 2-pack	0816-85-759



<b>Arm Measurement Guide</b>	
H	Acromion
G	Axilla
F	Midpoint of E and G
E	Elbow crease
D	Midpoint of C and E
C1	5cm above C
C	Wrist
B	Base of thumb webbing
A	Base of little finger
Z	Distal interphalangeal joint
ℓ = length	
c = circumference	