



TILT TABLE



350 lb capacity

Motorized

Item No.

TG2724

FEATURES:

- 350 lb Weight Capacity
- UL Recognized Components
- Casters - 5" Swivel with Locking Brakes
- Removable Footboard
- Safety Belts
- Angle Indicator
- Accommodates Patient Lifts
- 2" Upholstered Top
- Steel Frame

| | |
|-------------|---|
| Width: | 28" |
| Length: | 78" |
| Height: | 33" |
| Tilt Range: | 0 to 90 degrees |
| Base: | 3" Tubular Steel |
| Finish: | Pearl Grey Standard; Porcelain Optional. |
| Upholstery: | Imperial Blue, Regimental Blue, Wedgewood, Blue Ridge, Dove, Graphite, Black, Burgundy, Mauve, Forest, Grotto, Sand, Camel. |
| Electrical: | 120V, 1.6A, 192W, hospital grade plug. |
| Warranty: | Ten (10) Year Limited Warranty. |

*Nominal Dimensions Stated