



TILT TABLE

 350 lb capacity

Motorized, Wide Footboard

Item No.

TG2702

FEATURES:

- 350 lb Weight Capacity
- UL Recognized Components
- Wide Footboard for Patient Comfort
- Casters - 5" Swivel with Locking Brakes
- Safety Belts
- Angle Indicator
- Accommodates Patient Lifts
- 2" Upholstered Top
- Steel Frame

Width: 28"
Length: 78"
Height: 33"
Tilt Range: 0 to 90 degrees
Base: 3" Tubular Steel
Finish: Pearl Grey Standard; Porcelain Optional.
Upholstery: Imperial Blue, Regimental Blue, Wedgewood, Dove, Graphite, Black, Burgundy, Mauve, Forest, Grotto, Blue Ridge, Sand, Camel.
Electrical: 1.6A, 120V, 192W, hospital grade plug.
Warranty: Ten (10) Year Limited Warranty.

*Nominal Dimensions Stated