

INSTRUCTIONS




NEUTRAL POSITION WRIST/HAND/FINGER SPLINT

- Apply Splint
 - Adjust straps and position pads for a comfortable fit
- Adjust Strap Length
 - Forearm Straps- Remove strap · cut non hook end to length · re-apply
 - Hand Straps- Cut in front of stitch line
- Plastic can be modified with a heat gun- Remove from patient · remove liner · avoid heating straps
 - Using extreme caution, gently heat plastic until soft · reshape · allow splint to cool completely before reapplying liner
- Cleaning
 - Plastic- Use damp cloth
 - Liner- Remove liner · hand wash · air dry completely

Remove splint from patient to make adjustments

Avoid storing in extreme heat or cold.

Single Patient Use. This product is to be used under the supervision of a medical professional



Authorized Representative: Medica Surgical Innovations Limited, Moorgate Street, Blackburn BB2 4PB, Lancs UK

INSTRUCTIONS



NEUTRAL POSITION WRIST/HAND/FINGER SPLINT

- Apply Splint
 - Adjust straps and position pads for a comfortable fit
- Adjust Strap Length
 - Forearm Straps- Remove strap · cut non hook end to length · re-apply
 - Hand Straps- Cut in front of stitch line
- Plastic can be modified with a heat gun- Remove from patient · remove liner · avoid heating straps
 - Using extreme caution, gently heat plastic until soft · reshape · allow splint to cool completely before reapplying liner
- Cleaning
 - Plastic- Use damp cloth
 - Liner- Remove liner · hand wash · air dry completely

Remove splint from patient to make adjustments

Avoid storing in extreme heat or cold.

Single Patient Use. This product is to be used under the supervision of a medical professional



Authorized Representative: Medica Surgical Innovations Limited, Moorgate Street, Blackburn BB2 4PB, Lancs UK