

ROLYANFIT WRIST AND THUMB SPICA

Product Features:

- Durable lightweight fabric is lycra lined for breathability and patient comfort.
- Removable aluminum stays provide support for palmar surface of wrist and extensor surface of thumb. Stays can be adjusted for proper angulation for customized fit.
- D-ring closure straps for easy application.
- 8" length.
- Latex free.

Indications:

- Gamekeeper's thumb, DeQuervain's tendonitis, sprains, post-cast healing, selected scaphoid injuries and carpal tunnel syndrome.



Left	Right	Size	Wrist Circ.
9298-89	9298-90	XSmall	5¼"-6¼"
9272-24-02	9272-24-01	Small	6¼"-7¼"
9272-25-02	9272-25-01	Medium	7¼"-8"
9272-26-02	9272-26-01	Large	8"-9"
9272-27-02	9272-27-01	XLarge	9"-10"