



High Percentage Body Unweighting System

GENERAL DESCRIPTION

GlideTrak™ is a highly versatile and affordable Body Unweighting System for pain free exercise and gait training. It features a unique Pelvic Suspension Design to lift the body and give complete support between steps.

FEATURES AND BENEFITS

High Percentage Body Unweighting

- Innovative design allows partial to full Body Unweighting capability
- Simple strap or treadmill incline adjustments control desired amount of unweighting
- Meets doctors orders for:
Non Weight Bearing, Touch Down Weight Bearing, Partial (50% or less) Weight Bearing or Weight Bearing as Tolerated
- Challenges patients from beginning gait training to full exercise or performance training

Designed to Fit Over Any Treadmill

- GlideTrak™ converts YOUR treadmill to a High Percentage Body Unweighting System
- Standard safety harness for patient confidence
- Frame rolls away on locking wheels for over the ground training or storage
- Easily assembled or taken apart in minutes
- Remove saddle any time for regular treadmill use

Hands-On Access to Patients

- Unobstructed views and hands-on access to patients during treatment
- Facilitate gait corrections and maximize hip extension
- Integrate trunk and arm movements
- Provide assistance or resistance



Patented Ergonomic Pelvic Suspension

- Securely cradles pelvis and lifts body weight
- No pressure on sensitive perineal area
- Provides 100% body support between each step
- Lift one or both legs any time to rest, ease pain, or correct gait and let the saddle completely support body weight
- Train with a single leg or gradually integrate the involved limb (**Differential Training**)
- Easy Fit Seat quickly adjusts to each patient





High Percentage Body Unweighting System



FEATURES AND BENEFITS

Safe Exercise and Postural Control

Exercise safely and securely on the GlideTrak™ with its standard shoulder harness

The optional **Balance Harness Upgrade** improves security, balance and postural control. It allows safe regular treadmill use and over the ground training for geriatric and other balance challenged patients. For use with or without the GlideTrak™ saddle assembly.



Neurological Treatment Benefits

- Offers the benefits of Body Unweighting, supported by numerous studies
- GlideTrak™ allows patients to sustain gait corrections at *faster than normal walking speeds*
- Faster gait speeds are the key that drives brain neuroplasticity resulting in enhanced gait, motor control and neurological recovery for many patients

Enhanced Sports and Orthopedic Recoveries

- Earlier pain free interventions
- Safe training through injuries
- Maintenance of peak fitness levels during recovery

Great for Use in Clinics and Homes

- Practical and affordable for home programs and lifelong exercise
- For more information visit: www.glidetrak.com

Packaging Details

| GlideTrak™ BODY UNWEIGHTED SYSTEM | |
|---|--|
| Standard Equipment | Easy Fit Seat, Ratchet Adjustors, Locking Wheel Kit, Standard Safety Harness |
| Assembly | Quickly Snaps Together or Apart |
| Optional Balance Harness Upgrade | Recommended for Postural & Balance Control, Comfort & Additional Security |
| Max Weight | 350 Pounds |
| GlideTrak™ Dimensions | 89" H X 42.5" W X 52" L |
| Shipping Dimensions | (Box: 44" X 37" X 10") |
| Weight | 98 Pounds |

Applications

- Orthopedic Injuries and Fractures
- Neurological Conditions
- Geriatric Balance & Posture Training
- Sports Injury Rehabilitation
- Pre and Post Operative Conditioning
- Joint Replacements
- Joint pain and Arthritis
- Sports Performance Enhancement
- Fitness & Cross Training
- Prosthetic & Single Leg Training
- Weight Loss Intervention